# InertSustain / Inertsil シリーズ

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Packed Column for High-Performance Liquid Chromatography  
**InertSustain / Inertsil Series**  
**Instruction Manual**

## 1 Introduction

Thank you very much for purchasing this product.

Please read this instruction manual carefully and use the product properly to ensure the best performance.

## 2 Handling Precautions

### 2-1 Precautions

- Do not drop or hit the column. Strong impacts may cause column degradation.
- The columns provide high pressure resistance due to using a high-pressure slurry method. However, in order to ensure long-term stable use, we recommend using the column at lower pressure than the values shown in the table below.

Particle diameter	Column size	Recommended pressure
1.7 μm	All sizes	18000 psi or less
1.7 μm	UHPLC PEEK Column	14500 psi or less
1.9 μm, 2 μm	All sizes	11600 psi or less
3 μm HP	All sizes	7200 psi or less
3.5 μm (Hybrid-C18)	All sizes	7200 psi or less
5 μm (Hybrid-C18)	All sizes	5800 psi or less
3 ~ 10 μm	I.D. 0.05 to 50 mm	2900 psi or less
5, 10 μm	Inner diameter 100 mm	1450 psi or less


- Be careful of sudden pressure fluctuations.
  - When removing the column, wait until the pressure gauge indicates 0.
  - Please note that slow operation of the sample injection valve will cause sudden pressure changes at the column inlet.
- To the extent possible, dissolve samples in solvents with the same composition as eluents. (initial solvent for gradients). Injecting a large amount of samples dissolved in solvents with higher solubility than eluents may reduce the separating capacity and cause the samples to precipitate at the inlets of the columns.
- In the case of replacing the columns with eluents that do not mix with each other, (e.g., changing from a normal-phase type to reversed-phase type), please run 2-propanol or ethanol as intermediate solvents between them. Please run these solvents at a capacity at least 10 times larger than the column volumetric capacity. Running alcohol-based solvents tends to indicate high pressure, therefore, adjust the flow rate, paying attention to the upper pressure limit of the column.
- If a peak that elutes quickly is tailing, it may be caused by dead volume. Check if the connection piping at the column joint is inserted all the way to the end. For pipes connecting to injectors and detectors, please select pipes that are suitable for the inner diameter of the columns in use or correspond to analytical systems. Particularly, in case of analyzing at low flow rates using semi-micro columns, etc., the pipe specification may affect the analysis result more significantly.
- Pressure increase and peak cracking may be caused by clogging or contamination of the filter at the column inlet.
  - Filter the eluent through a membrane filter of 0.45 μm or less before use.
  - Filter the sample solution with a syringe filter of 0.45 μm or less before injecting.
  - Filter clogging can be prevented by using Guard Column for UHPLC or Cartridge Guard Column E.
- Ensure that the column is well equilibrated with the eluent before use. (In case of HILIC mode or reversed-phase mode using ion-pair reagents, it may take longer to stabilize.)
- For Inertsil Amide, do not use solvents with less than 50% organic content.
- Depending on LC/MS analysis conditions and types of equipment, precipitates may be seen at the column outlet. Check the effect of precipitation before use.
- Refer to the table below for the applicable pH range and upper temperature limit of the column. (Regarding not listed products in the table, refer to the GL Sciences website for the applicable pH range and upper temperature limit.)
- Please note that exceeding the upper temperature limit lead to premature deterioration.

Product name	Applicable pH range (Normal 20 to 40 °C)	Upper temperature limit	
Inertsil Hybrid-C18	1-12 ※1	80°C (pH 1-10)	60°C (pH 1-12)
InertSustain C18, Swift C18, Bio C18, AQ-C18	1-10 ※2,3,4	60°C (pH 1-7)	50°C (pH 1-10)
InertSustain C8, Swift C8, Phenylhexyl	1-10 ※2,3,4	60°C (pH 2-7)	50°C (pH 1-9)
InertSustain AX-C18	1-9 ※2,3,4	60°C (pH 2-7)	50°C (pH 1-8)
InertSustain Amide	2-8.5 ※3	60°C (pH 2-7)	50°C (pH 2-8.5)
InertSustain Phenyl, NH <sub>2</sub> , Cyano, PFP	2-7.5 ※3	60°C (pH 2-7)	50°C (pH 2-7.5)
Inertsil Series	2-7.5 ※3	60°C (pH 2-7)	50°C (pH 2-7.5)
InertSustain C30	1-7.5 ※3	60°C (pH 1-7)	50°C (pH 1-7.5)

- ※1 Column life varies greatly depending on conditions such as applicable pH, temperature, and eluent compositions. To ensure stable and long-term use, we recommend lowering the column temperature, using low concentrations of buffer salts and additives, and analyzing with eluents containing organic solvents. For analysis at pH 1-2, we recommend using TFA, formic acid, acetic acid, phosphoric acid, etc. For analysis at pH 9-12 within the applicable pH range, we recommend the use of inorganic or organic buffer solutions in 1 - 20 mM range, and avoid using sodium hydroxide as mobile phase.
- ※2 Column life varies greatly depending on conditions such as applicable pH, temperature, and eluent compositions. To ensure stable and long-term use, we recommend lowering the column temperature, using low concentrations of buffer salts and additives, and analyzing with eluents containing organic solvents. For analysis at pH 1-2, we recommend using TFA, formic acid, acetic acid, phosphoric acid, etc. For analysis at pH 9-10 within the applicable pH range, we recommend the use of organic buffer solutions of such as triethylamine with a concentration of around 5 mM. For analysis using only buffer solutions that do not contain organic solvents, use them within pH 2-8.
- ※3 In order to prevent premature deterioration, do not exceed the pH of the eluent described above.
- ※4 In case of using the column at pH 1-2 or 9-10, we recommend analyzing at low temperature or using eluents with organic solvents such as methanol.

## 2-2 UHPLC PEEK / PEEK Column Handling Precautions

- All wetted parts of the UHPLC PEEK / PEEK column are made of resin. The column port in particular has a very delicate structure. Therefore, if the column is handled in the same manner as a stainless steel column, the wetted parts may be damaged.
- Since it can be connected with less torque than a stainless steel column, when connecting to equipment piping, etc., make sure that there is no leakage of fluid, and tighten with sufficient torque (about 0.8 N·m) to prevent leakage.
- Do not use a male nut with a deformed tip that is used for prolonged periods of time or repeatedly. Since the column port may be damaged, we recommend using new resin male nuts when making connections. [Recommended fitting: 2-piece type resin male nut.
- Running THF or chloroform for prolonged periods of time may lead to premature degradation.
- When installing or removing the column, please hold the male nut and the column end nut portion, and proceed.

## 2-3 Capillary EX / Capillary EX-Nano Column Handling Precautions

- The column joint is a Parker type (UP type) with 1/16" male nut connection and 1/32" connection.
- Never bend the capillary EX column.
- The capillary EX-nano column is somewhat bendable, however, do not bend it too much as this can degrade the column.
- If the piping is not connected correctly for the capillary column, the column performance will be negatively impacted due to dead volume. The connection kit (sold separately) simplifies the connection.

## 2-4 Preparative HPLC Column Handling Precautions

- Since the flow rate of the mobile phase is considerably higher than that of the normal analytical column, use a flow pipe with an inner diameter of 0.8 mm or 1.0 mm according to the application. (Can be used with general-purpose piping up to 14 mm inner diameter.)
- Higher flow rates can cause the column to be clogged due to blockages caused by valve switching during sample injection. This can lead to column degradation. We recommend attaching a bypass to protect the column.
- Refer to the following table for the optimum flow rate.

Column inner diameter	Optimum flow rate
7.6 mm	2 ~ 4 mL / min
8.0 mm	2 ~ 4 mL / min
10 mm	3 ~ 5 mL / min
14 mm	5 ~ 10 mL / min
20 mm	10 ~ 20 mL / min
30 mm	20 ~ 45 mL / min
50 mm	70 ~ 130 mL / min
100 mm	200 ~ 300 mL / min

### 3 Product Features

This product is manufactured under our strict standards at each production process for matrix silica gel, chemical modification, and column performance inspection after packing.

This ensures that the product can be used at a consistent quality levels with peace of mind.

### 4 Confirmation of Contents

- Check for any abnormalities in column appearance, packaging, etc.
- Check for any errors in packing material name, column size, etc.
- Check to see that a column performance report is included.  
The column performance report includes information such as the packing material lot No., column serial No., and column performance test eluent. Please keep it in a safe location.
- In principle, columns are shipped containing solvents described in the column performance report.

### 5 Options

Product name	Cat. No.
Connection Kit for Capillary EX Column	5020-01880
Connection Kit for Capillary EX-Nano Column	5020-01881

### 6 Specifications

Type	Main column	Eluent sealed at shipping
Normal-phase type	InertSustain NH2 Inertsil SIL-100A, CN-3, NH2, Diol, etc.	Hexane / ethanol Mixed eluent
Reversed-phase type	Inertsil Hybrid-C18 InertSustain C18, AQ-C18, AX-C18, Phenyl, PFP, Cyano, C30, etc. Inertsil ODS-4, ODS-3, C8-3, Ph-3, etc.	Acetonitrile / water Mixed eluent
HILIC type	InertSustain Amide Inertsil HILIC, Amide, etc.	Acetonitrile / water Mixed eluent
	InertSustain NH2 with 100% CH3CN, etc.	Acetonitrile 100%
Ion exchange type	Inertsil CX, AX, etc.	Methanol 100%

### 7 Washing and Storage

- When using reversed-phase columns, thoroughly wash with salt-free eluents if the eluents contain buffer salts or ion-pair reagents after use.
- When strong ionic components are analyzed with InertSustain AX-C18, wash it thoroughly with a methanol solution containing 50 mM to 100 mM ammonium acetate.
- Store the reversed-phase columns after replacement with an organic solvent such as acetonitrile or methanol.
- NH<sub>2</sub> column: When using reversed-phase eluents, wash the columns with water/acetonitrile in a 50/50 ratio, etc.
- When using normal-phase columns, wash with ethanol or 2-propanol.  
Generally speaking, running alcohol-based solvents tends to indicate high pressure. Therefore, adjust the flow rate, paying attention to the upper pressure limit of the column and reduce flow rate if necessary. Store after replacement with 100% hexane.
- After using the HILIC type columns with eluents containing buffer salts, run them with eluents containing high water concentration (50% or higher water) to remove hydrophilic substances, and store them after replacing with a high acetonitrile concentration (80% or higher acetonitrile).
- For Inertsil CX and Inertsil AX, they are shipped containing 100% methanol. Therefore, first replace them with pure water before running them with buffer solutions. After using the buffer solutions, first replace them with pure water and then store them with 100% methanol.
- When storing the columns, tightly plug them with the provided plug and store them in a cool place where there is little temperature change and low humidity.  
\* For UHPLC PEEK / PEEK columns, mount the plug in the same manner as when connecting to equipment, while being careful not to overtighten them.

InertSustain / Inertsil Series products are manufactured, inspected, packaged, and shipped under strict quality control, but in the unlikely event of a defect, please contact one of our branches, sales offices, or distributors.

However, we do not provide any warranty against damage, problems related to service life, or deterioration caused by use of the product in a manner not in compliance with this instruction manual.

The specification and appearance are subject to change without notice.

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