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1 Introduction

Thank you for purchasing the IM column. To provide optimum performance of the IM column hardware and to ensure its long-term usage, please read the following instructions before use. For the specification on the packing materials, please refer to the attached InertSustain and Inertsil care and use manual. Provided that a part of the specifications between the column hardware and the packing material is different, please give priority to instructions of the IM column hardware over those of the packing material

2 Safety Precautions

- When analyzing compounds which show a specific adsorption on SUS surface, please note that there is the possibility of an additional adsorption to the wetted parts of analytical instruments, such as flow path and autosampler needle. Therefore, please thoroughly rinse the instruments before you use them. Generally, it is known that the specific adsorption of analytes on SUS surface in conventional LC instruments can be reduced by applying the passivation. A typical conditioning procedure is shown as below.^{※1}

※1 Passivation of LC instruments

- Flush the LC instruments with 50% CH₃CN containing a 0.1% phosphoric acid aqueous solution at a half of the applied flow rate to the sample analysis for 24 hours or longer, while keeping the column temperature ambient (room temp.) (Note: since phosphoric acid is a non-volatile acid, please disconnect the flow tube from a detector such as MS where contamination is concerned).
- After washing, please replace mobile phase and connect the column and flow tube to the detector for measurements.

If the performance of IM column may not be adequate for separation despite the compliance to the above-mentioned specifications and cautions, the below conditioning can be firstly recommended.^{※2}

※2 Column conditioning

Flush the column with 50% CH₃CN containing a 0.1% formic acid aqueous solution at a room temp., using a half of the applied flow rate to the sample analysis for 10 hours or longer.

3 Product Features

- The IM column hardware is prepared with a steel-use-stainless (SUS) body with InertMask^{※3} coating. This product enables to reduce a specific adsorption of analytes particularly on column hardware surface and leads to improving sensitivity and peak shape. The color of the IM column hardware is different from that of a conventional SUS column, but it was already proved by the manufacturers that there is no effect on the column performance.
- ※3 Inert layer coating technology to the SUS surface.

4 Specification

To provide the optimum performance, following specifications are recommended.

IM Column

Parameter	Recommended range ^{※4}
Operating Pressure (1.9 μm, 2 μm)	≤ 80 MPa
Operating Pressure (3 μm)	≤ 50 MPa
pH	2.0-7.5

IM Column Guard^{※5}

Parameter	Recommended range ^{※4}
Operating Pressure	≤ 80 MPa
pH	2.0-7.5

※4 For the long-term usage, we recommend applying an eluent containing organic solvent and a low concentration buffer at a low measurement temperature (40°C ≦). After measurements, please flush 50% CH₃CN prepared with a 0.1% formic acid aqueous solution into this product at room temp. for longer than 1 hour, for instance, by running a LC pump at a flow rate of 0.10 mL/min.

※5 Use holder for Guard Columns for UHPLC (Cat. No.5020-08630) when using IM Column Guard

5 Storage

- Regarding the column storage, please refer to the attached InertSustain and Inertsil instruction manual.
- Provided that one conducts the analysis of strong anions with the InertSustain AX-C18 stationary phase, after the measurement we recommend flushing the IM column with 80% CH₃OH containing 50–100 mM ammonium acetate (or a higher methanol concentration) at room temp for 30–60 min, using the same flow rate as that of the sample analysis. Then, the IM column can be stored with 50% CH₃CN prepared with a 0.1% formic acid aqueous solution according to the procedure show in "Specification".
- The IM columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact us or your GL Sciences' supplier, who will ensure you complete satisfaction. We regret that we cannot guarantee the lifetime of columns, nor can we accept any claim when their performance has deteriorated due to no-compliance with the above operating manual or as a result of normal aging.