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## ISFET pH Electrode 0040N-10D Instruction Manual

Thank you for purchasing a HORIBA ISFET pH electrode. This product is less fragile and more easy to handle than existing glass pH electrodes because of semiconductor, and can be used with pH meters (D-50, D-70, D-200, F-50 and F-70 series). Read this manual carefully before use. This is the English translation of an original Japanese document.

### ■ Warranty and responsibility

The product (the converter excluding batteries, sensor, and accessories) shall be repaired or replaced free of charge should any malfunction occur during the course of normal use if returned within the warranty period (12 months from date of purchase). Return the converter with proof of purchase date.

This warranty does not cover damages caused by accidents, disasters, misuse, carelessness, and unauthorized disassemble and remodeling.

The sole and exclusive remedy of HORIBA Advanced Techno Co., Ltd. under this limited warranty shall be repair or replacement of the product within a commercially reasonable time from notification.

In any case, HORIBA Advanced Techno Co., Ltd. shall not be responsible for customer's passive damages, or damages which may arise from claims instituted by a third party against customers.

### ■ EU regulations

#### ● Conformable standards

This equipment conforms to the following standards:

CE	EMC: EN61326-1 Class B, Basic electromagnetic environment
	RoHS: EN IEC 63000 9. Monitoring and control instruments including industrial monitoring and control instruments

**Warning:** This product is not intended for use in industrial environments. In an industrial environment, electromagnetic environmental effects may cause the incorrect performance of the product in which case the user may be required to take adequate measures.

#### ● Information on disposal of electrical and electronic equipment and disposal of batteries and accumulators


The crossed out wheeled bin symbol with underbar shown on the product or accompanying documents indicates the product requires appropriate treatment, collection and recycle for waste electrical and electronic equipment (WEEE) under the Directive 2012/19/EU, and/or waste batteries and accumulators under the Directive 2006/66/EC in the European Union.

The symbol might be put with one of the chemical symbols below. In this case, it satisfies the requirements of the Directive 2006/66/EC for the object chemical.

This product should not be disposed of as unsorted household waste.

Your correct disposal of WEEE, waste batteries and accumulators will contribute to reducing wasteful consumption of natural resources, and protecting human health and the environment from potential negative effects caused by hazardous substance in products.

Contact your supplier for information on applicable disposal methods.

- **Authorised Representative in EU**  
HORIBA Europe GmbH  
Hans-Mess-Str.6, D-61440 Oberursel, Germany

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### ■ Safety precautions

<b>CAUTION</b>
<p> Handle carefully the sharp electrode tip to not get hurt.</p>
<p> If calibration standard solution gets on the skin, wash off it completely. In case that the solution gets into the eyes, continue to flush the eyes with running water for at least 15 minutes. And after washing the eyes thoroughly, get medical help immediately.</p>

### ■ Points of concern

- Only use the product including accessories for their intended purpose.
- The materials used of the electrode may dissolve in the sample depending on the specified sample or lengthy measurement. Do not eat the whole, or the part around where the electrode sensor touched, of measured samples.  
If measured sample gets on the skin, wash off it completely with soap and water.
- Do not disassemble any parts not specified in this manual. This may cause malfunction.
- Do not touch the sensing part and sensor terminal, or It may be damaged by static electricity.
- Store the product in dry conditions. Storage area should be kept away from high temperature, high humidity, and direct sunlight.
- Do not drop the electrode. And do not apply excessive force to the sensor and converter.
- Avoid continuous use in high temperature, acids or alkali sample of 0.01 mol/L or thicker. Failure to comply may cause decline life period of the sensor.
- Samples containing strong reducing agents including acidic thiourea may deteriorate the sensing part of the sensor.
- Make sure to use the electrode at a temperature within the specified usable range (sensor: 0°C to 60°C, converter: 0°C to 40°C).
- The sensor is a consumable. If it is damaged or deteriorated, replace a new sensor. The sensor cannot be repaired.
- White powder or solution may be generated at the liquid junction. This is leaked internal solution of the reference electrode, and not a symptom of trouble. Rinse the sensor with ion-exchange water before use.
- If internal solution of the reference electrode gets on the skin, wash off it completely. In case that the solution gets into the eyes, continue to flush the eyes with running water. And after washing the eyes thoroughly, get medical help immediately.

- Black spots visible near the sensing part of the 0141N result from the epoxy adhesive used in the assembly process and do not affect measurements.

### ■ Precautions in measurements

- In the measurement of liquid sample, if the sample has bubbles, the bubbles may adhere to the sensing part. Bubbles on the sensing part may affect a measurement precision.  
If necessary, swing the sensor gently in the solution or lift the sensor from the solution and immerse it again to remove bubbles near the sensing part and liquid junction.
- Measure in a place avoiding direct sunlight and not by the window. If measurement is performed by the window, the measurement value may differ 0.1 pH.
- Dirt on the sensing part or liquid junction will affect accuracy adversely. Keep the sensor clean referring to “Cleansing the sensing part” (page 4)”.

### ■ Battery management

- The battery life is about 500 hours. The attached batteries may be shorter in life than new ones.
- Misusing batteries may cause an internal solution leak or a damage of the batteries. Follow the instruction below.
  - CR2032 batteries are not rechargeable.
  - Remove the batteries from the converter before prolonged nonuse.
  - If the battery solution leaks, wipe away the spilled solution from the battery case before replacing the batteries.
- For child’s safety, keep batteries away from children. If your child does swallow a battery, take him or her to the emergency room immediately.
- Do not throw batteries into fire.


## 1 Specifications and Components

### 1.1 Packaged Contents


Name	Q'ty
Sensor (0141N)	1 pc
Converter	1 pc
Battery CR2032	2 pcs
3.33 mol/L KCl solution (model 300, polyethylene container)	1 pc
Instruction manual (this book)	1 booklet
Screw driver	1 pc

### 1.2 Components

#### ● Sensor (0141N)

#### ● Converter

### 1.3 Specifications


Electrode model	0040N-10D
Sensor model	0141N
Measurement range	pH 0 to 14
Usable temperature range	[Sensor] 0°C to 60°C, [Converter] 0°C to 40°C
Storage temperature range	0°C to 40°C
Electrode length	190 mm (excluding cable)
Sensing device	ISFET (Ion-sensitive field-effect transistor)
Internal electrode	Silver/ silver chloride electrode
Reference electrode internal solution	Carboxymethyl cellulose (CMC), glycerin, saturated KCl
Liquid junction material	Porous sintered polyethylene
Sensor material	Tantalum pentoxide, LCP, epoxy resin, polyethylene resin, gold, titanium
Height of sensing part	1 mm (from the tip of electrode)
Application	D-50, D-70, D-200, F-50, F-70 series
Power	6 V DC CR2032 ×2, Auto power OFF (1 h)
LED indicator	Meas. LED (green blink while the power is ON), Low Batt. LED (red blink when the battery voltage is low)
Cable length	1 m

## 2 Preparation


For the first use, prepare the ISFET pH electrode as follows:

### 2.1 Setting batteries


1. Unhook and draw the converter cover slowly.

2. Insert 2 button batteries (CR2032) inside the converter with the (+) sign up. (Check images).

3. Attach and hook the converter cover to the converter by making the projection (A) fit in the hole (a) on the hook.


### 2.2 Assembling sensor and converter


#### Note

- Do not touch the sensing part and sensor terminal, or it may be damaged by static electricity.
- When tightening the nut, no gap should be left between the nut and the sensor. If the nut is not tightened enough, water may enter the inside of the sensor and the sensor may be damaged.


1. Remove the end cap from the sensor, and the seal from the liquid junction. Make sure that the O-ring is attached to the sensor.

2. Turn the nut to the left to remove it from the converter.

3. Insert the sensor to the converter slowly checking that the O-ring is untwisted, and match the  $\Delta$  marks on the sensor and converter. Then turn the nut to the right to fix the sensor.

### 2.3 Connecting to pH meter


Before calibration or measurement, prepare the ISFET pH electrode as follows:

1. Insert the electrode connector into the connector port on the meter, after aligning with the pin. Do not insert the connector unless it is aligned properly with the port.
2. Press the electrode connector into the port on the meter, while turning the connector to the right.
3. Slide the connector cover over the electrode connector. Then, push the cover in straight until it comes in light contact with the meter case. Do not turn the cover.
4. Insert the temperature connector securely into the jack on the meter.

### 3 Zero adjustment of the ISFET

ISFET needs the zero adjustment procedures before calibration with standard solution.


1. Immerse the sensing part and liquid junction completely into 3.33 mol/L KCl solution (model 300) for 10 minutes or longer.


#### Note

For the first use, or use after a long interval (one week or longer), immerse the sensor tip in 3.33 mol/L KCl solution (model 300) for approximately 30 minutes before use.

2. After that, rise them with ion-exchange water and wipe them by patting the electrode dry with filter paper.


3. Press the power switch to turn ON the power. The Meas. LED blinks green.


#### Note

Automatic power OFF works in approximately 1 hour after power ON.

4. Keep the electrode in pH7 standard solution, with the sensing part and liquid junction completely immersed. If necessary, swing the sensor gently in the solution or lift the sensor from the solution and immerse it again to remove bubbles near the sensing part and liquid junction.

5. Adjust the zero point of the ISFET pH electrode by turn the zero adjustment volume with a slot driver, to keep pH meter readings within the range shown below.

In mV measurement mode:  $0 \text{ mV} \pm 10 \text{ mV}$   
In pH measurement mode:  $\pm 0.2 \text{ pH}$  for the pH value of pH 7 standard solution


#### Note

- Do not turn the zero adjust volume over the range of movement, or damage the volume. Zero adjustment volume has a 250-degree range of movement.
- Before zero adjustment in pH measurement mode, clear the existing calibration values. See the instruction manual of the pH meter for the details.

6. Calibrate the ISFET pH electrode according to the instruction manual of the pH meter you use.

### 4 Measurement


1. Wash the sensing part and liquid junction with ion-exchange water and dry them by patting the electrode dry with filter paper, prior to every immersion into standard/sample solution.

2. Perform the following handling specified by sample, and measure the sample according to the instruction manual of the pH meter you use.


#### ● For liquid samples

Keep the electrode in the solution, with the sensing part and liquid junction completely immersed. If necessary, swing the electrode gently in the solution or lift the sensor from the solution and immerse it again to remove bubbles near the sensing part and liquid junction.


#### ● For semisolid sample such as agar

Stab the electrode in the sample until the sensing part and liquid junction are buried fully.


#### ● For solid sample such as meat

3. Press the electrode tip into the sample slowly.

4. Lean the electrode slowly to contact the sensing part and liquid junction with the sample. Make sure to keep the sensing part close to the sample without a gap.

- **For droplet on the surface of solid**  
Perform in the same manner for solid sample.

**Note**

The materials used of the electrode may dissolve in the sample depending on the specified sample or lengthy measurement. Do not eat the whole, or the part around where the electrode sensor touched, of measured samples.

If measured sample gets on the skin, wash off it completely with soap and water.


**Recommended conditions for accurate measurement**

For accurate measurement, note the items below.

- Stir the sample gently with a stirrer during measurement or calibration.
- Standard solutions for calibration and measurement samples should be the same in temperature. If temperature difference between them is 10°C measured value may differ 0.5 pH.
- Rinse the electrode with the sample before measurement.
- Calibrate the electrode before every measurement.
- Before use, immerse the sensor into 3.33 mol/L KCl solution (model 300) for 30 minutes or longer.

**5 Storage**

1. Turn OFF the power switch.  
The Meas. LED turns OFF.

2. Wash the electrode well with ion-exchange water to remove sample completely and wipe it with filter paper, and store it in dry conditions.

**Note**

Avoid storing the electrode in a high temperature and humidity place. Keep the electrode dry and out of direct sunlight.

**6 Maintenance**

**6.1 Cleansing the sensing part**

Dirt on the sensing part may affect the electrode response, sensitivity, or a measurement precision. If the electrode is very dirty and cannot be washed clean using ion-exchange water, clean it as follows.

1. Clean the sensing part and liquid junction with cleansing agent matched to its particular dirt.

Type of dirt	Cleansing agent	Cleaning procedure
General, Oily	Diluted neutral cleansing agent	<ul style="list-style-type: none"> <li>● Wipe off with gauze soaked in cleansing agent.</li> <li>● Scrape with a tooth brush soaked in cleansing agent.</li> </ul>
General, Oily	Alkalescent cleansing agent	Immerse the electrode in the cleansing agent. Set the immersing time accordingly depending on the dirt condition. Follow the instruction manual of each cleansing agent company how to use it.
Organic substance such as protein	Cleansing agent (model 250)	Immerse the electrode in the cleansing agent. See the instruction manual for cleansing agent for more information about immersing time, etc.
Silver chloride	5% ammonia water	Immerse the electrode in the cleansing agent for 30 minutes.
Inorganic substance	Hydrochloric acid less than 0.01 mol/L	Immerse the electrode in the cleansing agent for an hour.

**Note**

- Do not use organic solvent (e.g. acetone, THF), acids or alkali of 0.01 mol/L or more. Using these chemicals may cause decline life period of the sensor.
- Do not use the cleansing agent (model 220) for glass electrode. It may deteriorate the sensing part.

2. Rinse the sensing part and liquid junction with ion-exchange water completely.


**Note**

Be sure to calibrate the electrode after cleansing the sensing part.

**6.2 Replacing the sensor**

The sensor is a consumable. If any of the following symptoms is found after responsive-part cleaning, replace the sensor in the procedures below.

- Sensitivity or response speeds do not recover.
  - The electrode can not be adjusted in the way mentioned in "3 Zero adjustment of the ISFET" (page 3).
  - The electrode can not be calibrated with standard solutions.
1. Turn OFF the power switch.
  2. Turn the nut to the left to remove it. And pull out the sensor from the converter.

3. Attach a new sensor referring to "2.2 Assembling sensor and converter" (page 2).

**Note**

Do not touch the sensor terminal. It may be damaged by static electricity.

**6.3 Replacing batteries**

The battery life is about 500 hours.

When battery voltage is low, the Low Batt. LED blinks red. (The LED is off during battery shutoff.) In this case, replace batteries as follows:

1. Pull out the sensor from the converter referring to "6.2 Replacing the sensor" (page 4).
2. Replace batteries with new ones referring to "2.1 Setting batteries" (page 2). If used batteries are set, remove it by using a ball-point pen or the like.
3. Attach the sensor referring to "2.2 Assembling sensor and converter" (page 2).

**Note**

- Replace two batteries (CR2032) with new ones at once.
- Do not touch the sensor terminal. It may be damaged by static electricity.
- Make sure to remove the sensor prior to the converter cover. Removing the cover with the sensor attached causes malfunction.