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
ISFET pH Sensor 0141N Instruction Manual

Thank you for purchasing a HORIBA ISFET pH electrode. ISFET pH sensor (0141N) is a component of ISFET pH electrode (0040N-10D). Read this manual with that of ISFET pH electrode carefully before use.

1 Cautions on use

■ Safety Precautions

 CAUTION	
	Handle carefully the sharp electrode tip to not get hurt.
	If calibration standard solution gets on the skin, wash off it completely. In case that the solution gets into the eyes, continue to flush the eyes with running water for at least 15 minutes. And after washing the eyes thoroughly, get medical help immediately.

■ Points of concern


- The materials used of the electrode may dissolve in the sample depending on the specified sample or lengthy measurement. Do not eat the whole, or the part around where the electrode sensor touched, of measured samples.
If measured sample gets on the skin, wash off it completely with soap and water.
- Do not touch the sensing part and sensor terminal, or it may be damaged by static electricity.
- Store the product in dry conditions. Storage area should be kept away from high temperature, high humidity, and direct sunlight.
- Do not wet the electrode connector and the sensor terminal.
- Avoid continuous use in high temperature, acids or alkali sample of 0.01 mol/L or thicker. Failure to comply may cause decline life period of the sensor.
- White powder or solution may be generated at the liquid junction. This is leaked internal solution of the reference electrode, and not a symptom of trouble. Rinse the sensor with ion-exchange water before use.
- If internal solution of the reference electrode gets on the skin, wash off it completely. In case that the solution gets into the eyes, continue to flush the eyes with running water. And after washing the eyes thoroughly, get medical help immediately.
- Black spots visible near the sensing part of the 0141N result from the epoxy adhesive used in the assembly process and do not affect measurements.

2 Specifications and Components

2.1 Packaged Contents

Name	Q'ty
Sensor (0141N)	1 pc
Instruction manual (this book)	1 booklet

2.2 Components

2.3 Specifications

Sensor model	0141N
Storage temperature range	0°C to 40°C
Length	100 mm
Sensing device	ISFET (Ion-sensitive field-effect transistor)
Internal electrode	Silver/ silver chloride electrode
Reference electrode internal solution	Carboxymethyl cellulose (CMC), glycerin, saturated KCl
Liquid junction material	Porous sintered polyethylene
Sensor material	Tantalum pentoxide, LCP, epoxy resin, polyethylene resin
Height of sensing part	1 mm (from the tip-)
Applicable electrode	0040N-10D


3 Preparation

3.1 Assembling sensor and converter


Note

- Do not touch the sensing part and sensor terminal, or it may be damaged by static electricity.
- When tightening the nut, no gap should be left between the nut and the sensor. If the nut is not tightened enough, water may enter the inside of the sensor and the sensor may be damaged.

1. Remove the end cap from the sensor, and the seal from the liquid junction. Make sure that the O-ring is attached to the sensor.

2. Turn the nut to the left to remove it from the converter.

3. Insert the sensor to the converter slowly checking that the O-ring is untwisted, and match the Δ marks on the sensor and converter. Then turn the nut to the right to fix the sensor.