

NDK-2A-F (220V)  
SV=50°C 到達15分後測定


測定

Sp1	51.4 °C
Sp2	51.7 °C
Sp3	51.6 °C
Sp4	51.9 °C
Sp5	51.6 °C
Sp6	51.7 °C
Sp7	51.9 °C
Sp8	51.7 °C
Sp9	51.8 °C

パラメータ

放射率	0.92
反射温度	0 °C

2020/06/11 16:37:41



FLIR0269.jpg

FLIR E4

63983277

2020/06/11 16:37:41



FLIR0269.jpg

FLIR E4

63983277