


■白硬注射器 ロック・メタル先



単位(mm)																	
	容量	規格	A	B	C	D	E	ID	F	H	J	K	L	M	N	OD	Z
1-6187-01	1ml	－	72	<u>55.5</u> 57.0	<u>1.5</u> 1.7	<u>3.1</u> 3.6	<u>9.6</u> 10.2	6.5	<u>16.0</u> 17.0	<u>14.5</u> 15.5	<u>64.0</u> 66.0	<u>9.0</u> 9.7	<u>2.8</u> 3.8	<u>11.0</u> 12.0	<u>5.8</u> 6.9	6.5	≒85
1-6187-02	2ml	－	73	<u>56.0</u> 57.5	<u>1.2</u> 1.5	<u>3.1</u> 3.6	<u>12.1</u> 12.7	9.5	<u>17.0</u> 19.0	<u>15.0</u> 17.0	<u>66.0</u> 68.0	<u>11.0</u> 11.7	<u>3.0</u> 4.0	<u>13.8</u> 14.8	<u>7.4</u> 8.5	9.5	≒87
1-6187-03	3ml	－	83	<u>65.5</u> 67.0	<u>1.2</u> 1.4	<u>3.1</u> 3.6	<u>12.4</u> 13.0	10.1	<u>18.5</u> 19.5	<u>16.0</u> 17.0	<u>76.0</u> 78.0	<u>11.5</u> 13.0	<u>3.0</u> 4.0	<u>14.0</u> 15.0	<u>7.6</u> 8.8	10.1	≒96
1-6187-04	5ml	－	93	<u>74.0</u> 75.5	<u>1.3</u> 1.5	<u>3.2</u> 3.7	<u>15.1</u> 15.7	12.4	<u>21.5</u> 23.0	<u>20.0</u> 21.0	<u>85.0</u> 87.0	<u>12.0</u> 13.5	<u>3.2</u> 4.2	<u>16.7</u> 17.9	<u>9.4</u> 11.0	12.4	≒105
1-6187-05	10ml	－	115	<u>95.0</u> 97.0	<u>1.5</u> 1.7	<u>3.3</u> 4.0	<u>18.0</u> 18.6	15.0	<u>25.5</u> 27.0	<u>23.5</u> 24.5	<u>108.0</u> 110.0	<u>14.5</u> 16.0	<u>3.8</u> 4.9	<u>19.0</u> 20.5	<u>12.0</u> 13.8	15.0	≒130
1-6187-06	20ml	－	134	<u>114.5</u> 116.5	<u>1.8</u> 2.0	<u>4.3</u> 5.1	<u>22.6</u> 23.3	18.9	<u>32.0</u> 32.8	<u>29.0</u> 29.8	<u>131.0</u> 133.0	<u>17.5</u> 19.0	<u>4.3</u> 5.8	<u>24.0</u> 25.5	<u>15.0</u> 16.8	18.9	≒154
1-6187-07	30ml	－	140	<u>119.0</u> 121.5	<u>2.1</u> 2.3	<u>4.7</u> 5.5	<u>26.9</u> 27.7	22.6	<u>36.5</u> 38.5	<u>33.0</u> 35.0	<u>139.0</u> 142.0	<u>21.5</u> 23.5	<u>4.9</u> 6.1	<u>27.3</u> 28.8	<u>17.7</u> 19.5	22.6	≒163
1-6187-08	50ml	－	172	<u>152.0</u> 154.5	<u>2.3</u> 2.6	<u>5.5</u> 6.5	<u>30.0</u> 31.0	25.3	<u>41.5</u> 43.0	<u>39.0</u> 41.0	<u>173.5</u> 176.5	<u>23.0</u> 25.0	<u>5.9</u> 7.5	<u>31.5</u> 33.0	<u>20.8</u> 22.6	25.3	≒198
1-6187-09	100ml	－	182	<u>157.5</u> 161.0	<u>2.5</u> 2.7	<u>6.0</u> 7.0	<u>40.5</u> 41.9	35.7	<u>54.5</u> 56.5	<u>50.0</u> 52.5	<u>187.5</u> 190.5	<u>35.5</u> 38.1	<u>7.1</u> 9.2	<u>41.0</u> 43.0	<u>30.0</u> 32.5	35.7	≒215

受注生産(返品不可)

1-6187-10	200ml	-	239	<u>218.0</u> 220.0	<u>2.7</u> 2.9	<u>7.0</u> 8.0	<u>47.3</u> 48.7	42.4	<u>61.9</u> 63.0	<u>58.8</u> 60.2	<u>255.2</u> 256.6	<u>31.3</u> 32.7	<u>12.6</u> 13.0	<u>53.5</u> 54.6	<u>27.4</u> 28.8	42.4	≒279
	300ml	-	268	<u>242.0</u> 247.0	<u>3.0</u> 3.5	<u>8.0</u> 8.5	<u>52.0</u> 57.0	49.0	<u>66.0</u> 71.0	<u>64.5</u> 69.5	<u>277.0</u> 282.0	<u>32.0</u> 37.0	<u>13.0</u> 13.5	<u>56.0</u> 61.0	-	49.0	≒303
	500ml	-	286	<u>260.0</u> 265.0	<u>4.0</u> 4.5	<u>9.0</u> 9.5	<u>67.0</u> 72.0	62.0	<u>82.0</u> 87.0	<u>80.5</u> 85.5	<u>307.0</u> 312.0	<u>44.0</u> 49.0	<u>14.0</u> 14.5	<u>66.0</u> 71.0	-	62.0	≒333
	1000ml	-	298	<u>272.0</u> 277.0	<u>6.0</u> 6.5	<u>9.5</u> 10.0	<u>97.0</u> 102.0	88.0	<u>111.0</u> 116.0	<u>109.5</u> 114.5	<u>313.0</u> 318.0	<u>37.0</u> 42.0	<u>18.0</u> 18.5	<u>107.0</u> 112.0	-	88.0	≒338

※ガラス製の為、上記寸法値と異なる場合がございます。