

		SOCOREX アクキュラマイクロ855(AcuraMicro855)				SOCOREX アクキュラマイクロ855(AcuraMicro855)				
チップ	容量	商品コード	1-9349-01	1-9349-02	1-9349-03	1-9349-04	1-9350-01	1-9350-02	1-9350-03	1-9350-04
		型番	855.08.010	855.08.050	855.08.200	855.08.350	855.12.010	855.12.050	855.12.200	855.12.350
		分注容量	0.5-10 μ l	5-50 μ l	20-200 μ l	100-1000 μ l	0.5-10 μ l	5-50 μ l	20-200 μ l	100-1000 μ l
SOCOREX										
クオリテック	10 μ	1-5300-11(3090.0010B)	●	×	×	×	●	×	×	×
	10 μ	1-5300-12(3020.0010B)	×	×	×	×	×	×	×	×
	200 μ	1-5300-13(3070.0200B)	×	●	●	×	×	●	●	×
	350 μ	1-5300-04(308.0350B)	×	●	●	●	×	●	●	●
	1000 μ	1-5300-15(3070.1000B)	×	×	×	×	×	×	×	×
アズワン										
ピオラマクラチップ	10 μ	3-6504-01(V-10B) 3-6629-01(V-10R) 3-6632-01(V-10RS)	●	×	×	×	●	×	×	×
	200 μ	3-6504-02(V-200B) 3-6629-02(V-200R) 3-6632-02(V-200RS)	×	●	●	×	×	●	●	×
	1000 μ	3-6504-03(V-1000B) 3-6629-03(V-1000R) 3-6632-03(V-1000RS)	×	×	×	×	×	×	×	×
ピオラマクラチップ (フィルター付)	10 μ	3-6625-01(V-10F) 3-6509-01(V-10FS)	●	×	×	×	●	×	×	×
	200 μ	3-6625-02(V-200F) 3-6509-02(V-200FS)	×	×	●	▲	×	×	●	▲
	1000 μ	3-6625-03(V-1000F) 3-6509-03(V-1000FS)	×	×	×	×	×	×	×	×
アイビス										
アイビスピペットチップ	0.5-20 μ	2-3006-01(I-201C) 2-4064-04(I-296C)	●	×	×	×	●	×	×	×
		2-4584-01(IN122-201C) 2-3007-05(I-212C)	×	×	×	×	×	×	×	×
	1-200 μ	2-3006-02(I-503Y) 2-3009-01(I-596Y)	×	●	●	▲	×	●	●	▲
		2-4584-02(IN122-503Y) 2-3007-01(I-513Y)	×	×	×	×	×	×	×	×
	1-200 μ	2-3006-04(I-703Y) 2-3006-05(I-703C)	×	●	●	▲	×	●	●	▲
		2-3009-03(I-796Y) 2-4584-03(IN122-703Y) 2-3007-03(I-713C)	×	×	×	×	×	×	×	×
	1-250 μ	2-4064-02(I-801LC) 2-4064-05(I-896LC)	×	●	●	▲	×	●	●	▲
	1-300 μ	2-4584-06(IN122-801C)	×	×	●	▲	×	×	●	▲
100-1000 μ	2-3006-10(IN110-505C)	×	×	×	×	×	×	×	×	
100-1000 μ	2-3006-03(I-502B) 2-3009-02(I597-B)	×	×	×	×	×	×	×	×	
100-1200 μ	2-4584-07(IN122-502B) 2-3007-02(I-512B)	×	×	×	×	×	×	×	×	
アイビスフィルターチップ	0.1-10 μ	2-4064-03(IN110-804C) 2-4064-06(IN120-804C)	×	×	×	×	×	×	×	×
		2-4584-08(IN122-804C) 2-3007-06(IN112-814CH)	×	×	×	×	×	×	×	×
	1-20 μ	1-8548-01(IN110-805C) 1-8548-02(IN120-805C)	×	×	×	×	×	×	×	×
		2-4544-01(IN124-10S) 2-4545-01(IN126-10S)	●	×	×	×	●	×	×	×
	1-100 μ	2-4544-02(IN124-20S) 2-4545-02(IN126-20S)	×	▲	▲	▲	×	▲	▲	▲
		2-4544-03(IN124-100S) 2-4545-03(IN126-100S)	×	●	▲	▲	×	●	▲	▲
	1-200 μ	2-4544-05(IN124-200(L)S) 2-4545-05(IN126-200(L)S)	×	●	●	▲	×	●	●	▲
		2-4544-04(IN124-200S) 2-4545-04(IN126-200S)	×	×	●	▲	×	×	●	▲
100-1000 μ	2-4544-06(IN124-1KS1) 2-4545-06(IN126-1KS1)	×	×	×	×	×	×	×	×	
100-1200 μ	1-8548-03(IN124-1200S)	×	×	×	×	×	×	×	×	
エッペンドルフ										
ピペットチップ(epTip)	0.1-10 μ	2-4876-01(93401) 2-4874-01(93451)	●	×	×	×	●	×	×	×
		2-4875-01(93421) 2-4873-01(93461)	×	×	×	×	×	×	×	×
	0.1-20 μ	2-4878-01(93482)	×	×	×	×	×	×	×	×
		2-4876-02(93403) 2-4874-02(93453)	×	×	×	×	×	×	×	×
	0.5-20 μ	2-4875-02(93423) 2-4873-02(93463)	×	×	×	×	×	×	×	×
		2-4876-03(93404) 2-4878-02(93484)	×	●	●	▲	×	●	●	▲
	2-200 μ	2-4874-03(93454) 2-4875-03(93424)	×	●	●	▲	×	●	●	▲
		2-4873-03(93464) 2-4877-01(93474)	×	×	×	×	×	×	×	×
	20-300 μ	2-4876-04(93405) 2-4878-03(93485)	×	●	●	▲	×	●	●	▲
		2-4874-04(93455) 2-4875-04(93425)	×	×	×	×	×	×	×	×
50-1000 μ	2-4873-04(93465) 2-4876-05(93406)	×	×	×	×	×	×	×	×	
	2-4878-04(93486) 2-4874-05(93456)	×	×	×	×	×	×	×	×	
50-1250 μ	2-4875-05(93426) 2-4873-08(93466)	×	×	×	×	×	×	×	×	
	2-4877-02(93476) 2-4876-06(93407)	×	×	×	×	×	×	×	×	
ギルソン										
ダイヤモンドチップ	0.1-10 μ	1-6857-01(D10V) 1-6857-02(D10)	●	×	×	×	●	×	×	×
		1-6857-03(D10T)	×	×	×	×	×	×	×	×
		1-6857-04(DL10V) 1-6857-05(DL10)	×	×	▲	▲	×	×	▲	▲
	2-200 μ	1-6857-06(DL10T)	×	×	×	×	×	×	×	×
		1-6857-07(D200V) 1-6857-08(D200)	×	●	●	▲	×	●	●	▲
	100-1000 μ	1-6857-09(D200T) 1-6857-10(D1000V)	×	×	×	×	×	×	×	×
	1-6857-11(D1000) 1-6857-12(D1000T)	×	×	×	×	×	×	×	×	

▲はチップの最大容量まで