




































			R 1342	R 1345	R 2302	R 1381	R 1382	R 1385	R 1388	R 1389 (PTFEコート)	R 1311	R 1312	R 1313	R 1352	R 1355	R 1300	R 1302	R 1303	R 1375	R 1376	R 2311	
																						
			長さ (mm)	350	550	800	350	350	550	800	350	350	400	350	550	350	350	350	550	550	800	
			シャフト径 (mm)	8	8	13	8	8	10	10	8	8	10	8	8	8	10	8	8	10	13	
			羽根直径 (mm)	50	100	150	45	55	140	140	75	30	50	70	60/15	100/24	80	100	40	70	150	
			IKA	0000741000	0000741300	0000739000	0001296000	0001295900	0000477700	0000477800	0002343600	0002332900	0002333000	0002333100	0000756900	0001132700	0000513500	0002387900	0002746700	0000757700	0000757800	0000739500
			As-one	1-7327-01	4-2399-01	1-3152-03	1-3152-11	1-3152-07	1-3731-13	1-3731-14	61-0005-43	61-0005-36	1-7327-02	61-0005-37	1-3152-04	1-3731-15	1-3152-01	1-7327-03	61-0005-62	1-7327-04	1-3731-16	1-3731-17
NANOSTAR 7.5 digital	0025007173	4-2358-01	◎			◎		◎			◎									△		
MICROSTAR 7.5 digital	0025006557	3-8002-04	◎			◎		◎			◎									△		
MICROSTAR 7.5 control	0020103117	3-8002-11	◎			◎		◎			◎									△		
MICROSTAR 15 digital	0025006558	3-8002-05	◎			◎		◎			◎									△		
MICROSTAR 15 control	0020103118	3-8002-12	◎			◎		◎			◎									△		
MICROSTAR 30 digital	0025006559	3-8002-06	◎			◎		◎			◎									△		
MICROSTAR 30 control	0020103119	3-8002-13	◎			◎		◎			◎									△		
MINISTAR 20 digital	0025006564	65-1919-62	◎			◎		◎			◎									△		
MINISTAR 20 control	0025005587	63-1288-92	◎			◎		◎			◎									△		
MINISTAR 40 digital	0025006572	65-1919-63	◎			◎		◎			◎									△		
MINISTAR 40 control	0025005599	63-1288-93	◎			◎		◎			◎									△		
MINISTAR 80 digital	0025006580	65-1919-64	◎			◎		◎			◎									△		
MINISTAR 80 control	0025005612	63-1288-94	◎			◎		◎			◎									△		
EUROSTAR 20 digital	0004442004	1-9943-22	◎			◎		◎			◎									△		
EUROSTAR 40 digital	0004444004	1-9943-21	◎			◎		◎			◎									△		
EUROSTAR 60 digital	0004446004	1-3150-31	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 60 control	0004440004	1-3150-42	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 100 digital	0004238104	1-7326-25	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 100 control	0004028504	1-7326-36	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 200 digital	0003990004	1-7326-21	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 200 control	0003992004	1-7326-33	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 400 digital	002500753	63-5534-67	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎		◎
EUROSTAR 400 control	002500760	63-5534-68	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎		◎
EUROSTAR 20 high speed digital	0004028604	1-9944-21																				
EUROSTAR 20 high speed control	0025001466	63-1825-88																				
EUROSTAR 100 P4 control	0025005138	64-9092-02	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
EUROSTAR 200 P4 control	0004090004	1-7326-34	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
RW 20 digital	0025003596	1-9945-01	◎			◎		◎			◎		◎		◎		◎		◎	◎		
RW 28 digital	0005040004	1-6608-11	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎		
RW 47 digital	0004050009 (50Hz) 0004050010 (60Hz)	61-4948-44 61-4948-43		◎	◎			◎			◎		◎		◎		◎		◎	◎		◎

			R 1330	R 1331	R 1333	R 3000.1	R 3001.1	R 3003	R 3003.1	R 3003.2	R 3004	R 3004.1	R 3004.2	R 1001	R 1002	R 1401	R 1405	R 1402			
																					
			長さ (mm)	350	350	550	565	575	350	550	800	359	565	819	160	140	EUROSTAR 20 high speed digital/control本体付属の軸に取り付け				
			シャフト径 (mm)	8	8	10	10	10	8	10	10	8	10	10	4	4	ホモジナイザー-T60用R50機種シャフトに取り付け				
			羽根直径 (mm)	45	90	150	100	100	50	100	150	30	50	70	34	12	55	45	40		
			IKA	0002022300	0002022400	0002747400	0020001192	0020001195	0020001203	0020001204	0020001205	0020001206	0020001207	0020001208	0000527400	0000527500	0001242900	0001289800	0001243300		
			As-one	1-3152-09	1-7327-05	1-3731-18	63-1288-95	63-1288-96	63-1288-97	63-1288-98	63-1288-99	63-1289-01	63-1289-02	63-1289-03	-	-	1-9944-11	1-9944-12	1-9944-13		
NANOSTAR 7.5 digital	0025007173	4-2358-01	◎								△			◎	◎						
MICROSTAR 7.5 digital	0025006557	3-8002-04	◎								△			◎	◎						
MICROSTAR 7.5 control	0020103117	3-8002-11	◎								△			◎	◎						
MICROSTAR 15 digital	0025006558	3-8002-05	◎								△			◎	◎						
MICROSTAR 15 control	0020103118	3-8002-12	◎								△			◎	◎						
MICROSTAR 30 digital	0025006559	3-8002-06	◎								△			◎	◎						
MICROSTAR 30 control	0020103119	3-8002-13	◎								△			◎	◎						
MINISTAR 20 digital	0025006564	65-1919-62	◎								△			◎	◎						
MINISTAR 20 control	0025005587	63-1288-92	◎								△			◎	◎						
MINISTAR 40 digital	0025006572	65-1919-63	◎								△			◎	◎						
MINISTAR 40 control	0025005599	63-1288-93	◎								△			◎	◎						
MINISTAR 80 digital	0025006580	65-1919-64	◎								△			◎	◎						
MINISTAR 80 control	0025005612	63-1288-94	◎								△			◎	◎						
EUROSTAR 20 digital	0004442004	1-9943-22	◎			◎		◎			◎			◎	◎						
EUROSTAR 40 digital	0004444004	1-9943-21	◎			◎		◎			◎			◎	◎						
EUROSTAR 60 digital	0004446004	1-3150-31	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 60 control	0004440004	1-3150-42	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 100 digital	0004238104	1-7326-25	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 100 control	0004028504	1-7326-36	◎	◎		◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 200 digital	0003990004	1-7326-21	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 200 control	0003992004	1-7326-33	◎	◎	◎	◎		◎			◎		◎		◎		◎		◎	◎	
EUROSTAR 400 digital	002500																				