

【ピペットチップ - 適合メーカー一覧表】

容量	品番	メーカー						
		バイオヒット	エッペンドルフ	フィンピペット	ギルソン	ニチリョウ	レイニン	
マイクロチップ	10 $\mu$ l	521000/521100/522100/521105/521127	●	●	●	●	●	●
	10(20) $\mu$ l Extra Long	521050/521108/522108/521110/521112	●	●	●	●	●	●
	200 $\mu$ l	521014/521101/522101/521014-Y/521101-Y/522101-Y/ 521010/521010-Y/521106/521113	●	●	●	●	●	●
	250 $\mu$ l	521013/521107/522109	●	●	●	●	●	●
	300 $\mu$ l	521015/521102/522102	●	●	●	●	●	●
	1000 $\mu$ l	521016/521103/522103/521016-B/521103-B/522103-B/ 521020/521020-B/521017/521126/521115		●	●	●		
	1250 $\mu$ l	521040/521104/522104		●	●	●		
	5ml	521031		●				
		521032				●		
	10ml	521033		●		●		
	10ml Extra Long	521034		●		●		
	マキシペンス低滞留チップ	10 $\mu$ l	523100/524100/526100/521116/521121	●	●	●	●	●
10(20) $\mu$ l Extra Long		523109/524108/526108/521117/521122	●	●	●	●	●	●
200 $\mu$ l		523102/524101/526101/523101/521118/521123	●	●	●	●	●	●
250 $\mu$ l		523108/524105/526109	●	●	●	●	●	●
300 $\mu$ l		523103/524102/526102	●	●	●	●	●	●
1000 $\mu$ l		523105/524103/526103/523104/523106/521120/521125		●	●	●		
1250 $\mu$ l		523107/524104/526104		●	●	●		
フィルターチップ	10 $\mu$ l	527100/528100	●	●	●	●	●	●
	10(20) $\mu$ l Extra Long	527108/528108	●	●	●	●	●	●
	20 $\mu$ l	527101/528101	●	●	●	●	●	●
	50 $\mu$ l	527102/528102	●	●	●	●	●	●
	100 $\mu$ l	527103/528103	●	●	●	●	●	●
	200 $\mu$ l	527104/528104	●	●	●	●	●	●
	300 $\mu$ l	527105/528105	●	●	●	●	●	●
	1000 $\mu$ l	527106/528106		●	●	●	●	●
フィルターチップ・低滞留タイプ	10 $\mu$ l	529100	●	●	●	●	●	●
	10(20) $\mu$ l Extra Long	529108	●	●	●	●	●	●
	20 $\mu$ l	529101	●	●	●	●	●	●
	50 $\mu$ l	529102	●	●	●	●	●	●
	100 $\mu$ l	529103	●	●	●	●	●	●
	200 $\mu$ l	529104	●	●	●	●	●	●
	300 $\mu$ l	529105	●	●	●	●	●	●
	1000 $\mu$ l	529106		●	●	●	●	●