



商品名

アズフロン®ハーフネジチューブジョイント(オス)



品番	竹の子外径 ($C\phi$)	接続ネジ (R)	A mm	B mm	L mm
2-488-01	4	R1/8	14.0	11.0	31.0
2-488-02	6	R1/4	20.0	12.0	38.0
2-488-03	8	R1/4	20.0	12.0	38.0
2-488-04	10	R1/4	23.0	13.0	43.0
2-488-05	10	R3/8	23.0	13.0	43.0
2-488-06	12	R3/8	23.0	13.0	43.0

商品名	アズフロン®ハーフネジチューブジョイント(メス)
-----	--------------------------



品番	竹の子外径 (C φ)	接続ネジ (RC)	A mm	B mm	D mm	L mm
2-489-01	4	RC1/8	15.0	16.0	16.5	31.0
2-489-02	6	RC1/4	20.0	17.0	20.0	37.0
2-489-03	8	RC1/4	21.0	17.0	20.0	38.0
2-489-04	10	RC1/4	23.0	19.0	20.0	42.0
2-489-05	10	RC3/8	25.0	24.0	30.0	49.0
2-489-06	12	RC3/8	25.0	24.0	30.0	49.0