



コード	接続ネジ	サイズ (mm)					
		A	B1	B2	C	D	E
2-727-01	R 1/8×1/8	24	9	9	6	4	12
2-727-02	R 1/4×1/4	32	13	13	6	6	14
2-727-03	R 3/8×3/8	32	13	13	6	10	19
2-727-04	R 1/2×1/2	40	17	17	6	12	23
2-727-05	R 3/4×3/4	44	19	19	6	12	27
2-727-06	R 1/4×1/8	28	13	9	6	4	14
2-727-07	R 3/8×1/4	32	13	13	6	6	17
2-727-08	R 1/2×1/4	36	17	13	6	8	23
2-727-09	R 1/2×3/8	36	17	13	6	10	23

	UNIT	SCALE	DATE	材質	部品名
	mm	2:1			
承認	検図	製図		色	
				<b>アズワン株式会社</b>	
					図面番号



コード	接続径	サイズ (mm)			
		A	B	C	D
2-728-01	R 1/8×1/8	15	6	4	12
2-728-02	R 1/4×1/4	20	7	6	14
2-728-03	R 3/8×3/8	22.5	9.5	10	19
2-728-04	R 1/2×1/2	28.5	11.5	12	23
2-728-05	R 3/4×3/4	32.5	27	12	27
2-728-06	R 1/4×1/8	16	7	4	14
2-728-07	R 3/8×1/4	22.5	9.5	6	19
2-728-08	R 1/2×1/4	24.5	11.5	8	23
2-728-09	R 1/2×3/8	24.5	11.5	10	23

	UNIT	SCALE	DATE	材質	部品名
	mm	2:1			
承認	検図	製図			
				<b>アズワン株式会社</b>	
					図面番号



コード	接続径	サイズ (mm)						
		A	B	C	D	E	F	G
2-729-01	R 1/8	19	10	15	6	4	12	20
2-729-02	R 1/4	24.5	11.5	21.5	8.5	6	17	23
2-729-03	R 3/8	26.5	13.5	22.5	9.5	10	19	27
2-729-04	R 1/2	32	15	28.5	11.5	12	23	30
2-729-05	R 3/4	37	18	32.5	13.5	12	27	36

	UNIT	SCALE	DATE	材質	部品名
	mm	2 : 1			
承認	検図	製図			図面番号