

		SOCOREX 電動ピペット926XS (Acura electro pipettes926XS)										
チップ	容量	商品コード	3-5232-01	3-5232-02		3-5232-03	3-5232-04	3-5232-05	3-5232-06	3-5232-07		
		型番	926.0002PIU	926.0010PIU	926.0010PIU	926.0020PIU	926.0050PIU	926.0100PIU	926.0200PIU	926.1000PIU		
		商品コード	3-1609-01	3-1609-02		3-1609-03	3-1609-04	3-1609-05	3-1609-06	3-1609-07		
		型番	926.0002	926.0010	926.0010Y	926.0020	926.0050	926.0100	926.0200	926.1000		
		分注容量	0.1-2 μ l	0.5-10 μ l	1-10 μ l	2-20 μ l	5-50 μ l	10-100 μ l	20-200 μ l	100-1000 μ l		
SOCOREX												
クオリチップ	10 μ	1-5300-11(3090.0010B)	●	●	×	×	×	×	×	×		
	10 μ	1-5300-12(3020.0020B)	●	●	×	×	×	×	×	×		
	200 μ	1-5300-13(3070.0200B)	×	×	●	●	●	●	●	×		
	350 μ	1-5300-04(308.0350B)	×	×	×	×	●	●	●	×		
	1000 μ	1-5300-15(3070.1000B)	×	×	×	×	×	×	×	●		
アスワン												
ピオラモサクラチップ	10 μ	3-6504-01(V-10B) 3-6629-01(V-10R) 3-6632-01(V-10RS)	●	●	×	×	×	×	×	×		
	200 μ	3-6504-02(V-200B) 3-6629-02(V-200R) 3-6632-02(V-200RS)	×	×	●	●	●	●	●	×		
	1000 μ	3-6504-03(V-1000B) 3-6629-03(V-1000R) 3-6632-03(V-1000RS)	×	×	×	×	×	×	×	●		
ピオラモサクラチップ (フィルター付)	10 μ	3-6625-01(V-10F) 3-6509-01(V-10FS)	●	●	×	×	×	×	×	×		
	200 μ	3-6625-02(V-200F) 3-6509-02(V-200FS)	×	×	×	×	×	×	●	×		
	1000 μ	3-6625-03(V-1000F) 3-6509-03(V-1000FS)	×	×	×	×	×	×	×	●		
アイビス												
アイビスピペットチップ	0.5-20 μ	2-3006-01(I-201C) 2-4064-04(I-296C) 2-4584-01(IN122-201C) 2-3007-05(I-212C)	●	●	×	×	×	×	×	×		
	1-200 μ	2-3006-02(I-503Y) 2-3009-01(I-596Y) 2-4584-02(IN122-503Y) 2-3007-01(I-513Y)	×	×	●	●	●	●	●	×		
	1-200 μ	2-3006-04(I-703Y) 2-3006-05(I-703C) 2-3009-03(I-796Y) 2-4584-03(IN122-703Y) 2-3007-03(I-713C)	×	×	●	●	●	●	●	×		
	1-250 μ	2-4064-02(I-801LC) 2-4064-05(I-896LC) 2-4584-06(IN122-801C)	×	×	●	●	●	●	●	×		
	1-300 μ	2-3006-10(IN110-505C)	×	×	×	×	×	×	●	×		
	100-1000 μ	2-3006-03(I-502B) 2-3009-02(I597-B) 2-4584-07(IN122-502B) 2-3007-02(I-512B)	×	×	×	×	×	×	×	●		
	100-1000 μ	2-4064-03(IN110-804C) 2-4064-06(IN120-804C) 2-4584-08(IN122-804C) 2-3007-06(IN112-814CH)	×	×	×	×	×	×	×	●		
アイビスフィルターチップ	0.1-10 μ	1-8548-01(IN110-805C) 1-8548-02(IN120-805C)	×	×	×	×	×	×	×	●		
	1-20 μ	2-4544-01(IN124-10S) 2-4545-01(IN126-10S)	●	●	×	×	×	×	×	×		
	1-100 μ	2-4544-02(IN124-20S) 2-4545-02(IN126-20S)	×	×	●	●	▲	▲	▲	×		
	1-100 μ	2-4544-03(IN124-100S) 2-4545-03(IN126-100S)	×	×	×	×	●	●	▲	×		
	1-200 μ	2-4544-05(IN124-200(L)S) 2-4545-05(IN126-200(L)S)	×	×	×	×	●	●	●	×		
	1-200 μ	2-4544-04(IN124-200S) 2-4545-04(IN126-200S)	×	×	×	×	×	×	●	×		
	100-1000 μ	2-4544-06(IN124-1KS1) 2-4545-06(IN126-1KS1)	×	×	×	×	×	×	×	●		
100-1200 μ	1-8548-03(IN124-1200S)	×	×	×	×	×	×	×	●			
エッペンドルフ												
ピペットチップ(epTip)	0.1-10 μ	2-4876-01(93401) 2-4874-01(93451) 2-4875-01(93421) 2-4873-01(93461)	●	●	×	×	×	×	×	×		
	0.1-20 μ	2-4878-01(93482)										
	0.5-20 μ	2-4876-02(93403) 2-4874-02(93453) 2-4875-02(93423) 2-4873-02(93463)	●	●	×	×	×	×	×	×		
	2-200 μ	2-4876-03(93404) 2-4878-02(93484) 2-4874-03(93454) 2-4875-03(93424) 2-4873-03(93464) 2-4877-01(93474)	×	×	●	●	●	●	●	×		
	20-300 μ	2-4876-04(93405) 2-4878-03(93485) 2-4874-04(93455) 2-4875-04(93425) 2-4873-04(93465)	×	×	●	●	●	●	●	×		
	50-1000 μ	2-4876-05(93406) 2-4878-04(93486) 2-4874-05(93456) 2-4875-05(93426) 2-4873-08(93466) 2-4877-02(93476)	×	×	×	×	×	×	×	●		
	50-1250 μ	2-4876-06(93407) 2-4874-06(93457) 2-4875-06(93427) 2-4873-06(93467)	×	×	×	×	×	×	×	●		
	ギルソン											
	ダイヤモンドチップ	0.1-10 μ	1-6857-01(D10V) 1-6857-02(D10) 1-6857-03(D10T)	●	●	×	×	×	×	×	×	
		0.1-20 μ	1-6857-04(DL10V) 1-6857-05(DL10) 1-6857-06(DL10T)	×	×	×	×	×	×	×	×	
2-200 μ		1-6857-07(D200V) 1-6857-08(D200) 1-6857-09(D200T)	×	×	●	●	●	●	●	×		
100-1000 μ		1-6857-10(D1000V) 1-6857-11(D1000) 1-6857-12(D1000T)	×	×	×	×	×	×	×	×		

▲はチップの最大容量まで

		SOCOREX 電動ピペット936			
チップ	容量	商品コード	3-5232-08	3-5232-09	3-5232-10
		型番	936.02PIU	936.05PIU	936.10PIU
		商品コード	3-1609-08	3-1609-09	3-1609-10
		型番	936.02	936.05	936.10
		分注容量	0.2~2mL	0.5~5mL	1~10mL

SOCOREX

クオリチップ	2000 μ	1-9348-13(312.02)	●	×	×
	5000 μ	1-9348-11(312.05B)	×	●	×
	10000 μ	1-9348-12(312.10)	×	×	●

アイビス

アイビスピペットチップ	1000-5000 μ	2-3006-09(I-5000C)	×	●	×
-------------	-----------------	--------------------	---	---	---

エッペンドルフ

ピペットチップ(epTip)	250-2500 μ	2-4876-07(93408)	●	×	×
		2-4878-05(93488)			
		2-4874-07(93458)			
		2-4875-07(93428)			
	2-4873-07(93468)				
100-5000 μ	2-4876-08(93409)	×	●	×	
	2-4878-06(93489)				
	2-4874-08(93459)				
1000-10000 μ	2-4876-09(93400)	×	×	●※	
1000-10000 μ	2-4876-10(93410)	×	×	●※	

ギルソン

ダイヤモンドチップ	500-5000 μ	1-6857-13(D5000)	×	●※	×
-----------	----------------	------------------	---	----	---

※取付けられるがインジェクターが使用できない

作成2017年6月

チップ	容量	SOCOREX 電動ピベット956(Acura electro956)				SOCOREX 電動ピベット956(Acura electro956)				
		商品コード	3-5232-11	3-5232-12	3-5232-13	3-5232-14				
		型番	956.08.010PIU	956.08.050PIU	956.08.200PIU	956.08.350PIU	956.12.010PIU	956.12.050PIU	956.12.200PIU	956.12.350PIU
		商品コード	3-1609-11	3-1609-12	3-1609-13	3-1609-14				
		型番	956.08.010	956.08.050	956.08.200	956.08.350	956.12.010	956.12.050	956.12.200	956.12.350
		分注容量	0.5-10 μl	5-50 μl	20-200 μl	100-1000 μl	0.5-10 μl	5-50 μl	20-200 μl	100-1000 μl

SOCOREX

クオリチップ	10 μ	1-5300-11(3090.0010B)	●	×	×	×	●	×	×	×
	10 μ	1-5300-12(3020.0010B)	×	×	×	×	×	×	×	×
	200 μ	1-5300-13(3020.0200B)	×	●	●	×	×	●	●	×
	350 μ	1-5300-04(308.0350B)	×	●	●	×	×	●	●	●
	1000 μ	1-5300-15(3070.1000B)	×	×	×	×	×	×	×	×

アスワン

ピオラモサクラチップ	10 μ	3-6504-01(V-10B) 3-6629-01(V-10R) 3-6632-01(V-10RS)	●	×	×	×	●	×	×	×
	200 μ	3-6504-02(V-200B) 3-6629-02(V-200R) 3-6632-02(V-200RS)	×	●	●	×	×	●	●	×
	1000 μ	3-6504-03(V-1000B) 3-6629-03(V-1000R) 3-6632-03(V-1000RS)	×	×	×	×	×	×	×	×
ピオラモサクラチップ (フィルター付)	10 μ	3-6625-01(V-10F) 3-6509-01(V-10FS)	●	×	×	×	●	×	×	×
	200 μ	3-6625-02(V-200F) 3-6509-02(V-200FS)	×	×	●	▲	×	×	●	▲
	1000 μ	3-6625-03(V-1000F) 3-6509-03(V-1000FS)	×	×	×	×	×	×	×	×

アイビス

アイビスピベットチップ	0.5-20 μ	2-3006-01(I-201C) 2-4064-04(I-296C) 2-4584-01(IN122-201C) 2-3007-05(I-212C)	●	×	×	×	●	×	×	×
		2-3006-02(I-503Y) 2-3009-01(I-596Y) 2-4584-02(IN122-503Y) 2-3007-01(I-513Y)	×	●	●	▲	×	●	●	▲
	1-200 μ	2-3006-04(I-703Y) 2-3006-05(I-703C) 2-3009-03(I-796Y) 2-4584-03(IN122-703Y) 2-3007-03(I-713C)	×	●	●	▲	×	●	●	▲
		2-4064-02(I-801LC) 2-4064-05(I-896LC) 2-4584-06(IN122-801C)	×	●	●	▲	×	●	●	▲
	1-300 μ	2-3006-10(IN110-505C)	×	×	●	▲	×	×	●	▲
	100-1000 μ	2-3006-03(I-502B) 2-3009-02(I597-B) 2-4584-07(IN122-502B) 2-3007-02(I-512B)	×	×	×	×	×	×	×	×
	100-1000 μ	2-4064-03(IN110-804C) 2-4064-06(IN120-804C) 2-4584-08(IN122-804C) 2-3007-06(IN112-814CH)	×	×	×	×	×	×	×	×
		1-8548-01(IN110-805C) 1-8548-02(IN120-805C)	×	×	×	×	×	×	×	×
アイビスフィルターチップ	0.1-10 μ	2-4544-01(IN124-10S) 2-4545-01(IN126-10S)	●	×	×	×	●	×	×	×
		2-4544-02(IN124-20S) 2-4545-02(IN126-20S)	×	▲	▲	▲	×	▲	▲	▲
	1-100 μ	2-4544-03(IN124-100S) 2-4545-03(IN126-100S)	×	●	▲	▲	×	●	▲	▲
		2-4544-05(IN124-200(L)S) 2-4545-05(IN126-200(L)S)	×	●	●	▲	×	●	●	▲
	1-200 μ	2-4544-04(IN124-200S) 2-4545-04(IN126-200S)	×	×	●	▲	×	×	●	▲
		2-4544-06(IN124-1KS1) 2-4545-06(IN126-1KS1)	×	×	×	×	×	×	×	×
	100-1000 μ	1-8548-03(IN124-1200S)	×	×	×	×	×	×	×	×

エッペンドリフ

ピベットチップ(ep Tip)	0.1-10 μ	2-4876-01(93401) 2-4874-01(93451) 2-4875-01(93421) 2-4873-01(93461)	●	×	×	×	●	×	×	×
		2-4878-01(93482)								
	0.5-20 μ	2-4876-02(93403) 2-4874-02(93453) 2-4875-02(93423) 2-4873-02(93463)	×	×	×	×	×	×	×	×
		2-4876-03(93404) 2-4878-02(93484) 2-4874-03(93454) 2-4875-03(93424) 2-4873-03(93464) 2-4877-01(93474)	×	●	●	▲	×	●	●	▲
	20-300 μ	2-4876-04(93405) 2-4878-03(93485) 2-4874-04(93455) 2-4875-04(93425) 2-4873-04(93465)	×	●	●	▲	×	●	●	▲
		2-4876-05(93406) 2-4878-04(93486) 2-4874-05(93456) 2-4875-05(93426) 2-4873-08(93466) 2-4877-02(93476)	×	×	×	×	×	×	×	×
	50-1000 μ	2-4876-06(93407) 2-4874-06(93457) 2-4875-06(93427) 2-4873-06(93467)	×	×	×	×	×	×	×	×

キルソン

ダイヤモンドチップ	0.1-10 μ	1-6857-01(D10V) 1-6857-02(D10) 1-6857-03(D10T)	●	×	×	×	●	×	×	×
		1-6857-04(DL10V) 1-6857-05(DL10) 1-6857-06(DL10T)	×	×	▲	▲	×	×	▲	▲
	2-200 μ	1-6857-07(D200V) 1-6857-08(D200) 1-6857-09(D200T)	×	●	●	▲	×	●	●	▲
		1-6857-10(D1000V) 1-6857-11(D1000) 1-6857-12(D1000T)	×	×	×	×	×	×	×	×

▲はチップの最大容量まで