



ロングハンドルサンプラー サイズ詳細


5355-0120
5355-0250



5355-0300
5355-0450



5355-0180
5355-0280
5355-0600



品番	型番	部位	外径	長さ	重量 (kg)
3-9798-01	5355-0120	Rod①	Ø 20 mm	66 mm	0.2945
		Rod②	Ø 25 mm	60 mm	
		Joint①	Ø 25 mm	37,5 mm	
		Handle	Ø 37,5 mm	115,1 mm	
3-9798-02	5355-0180	Rod①	Ø 20 mm	55 mm	0.595
		Rod②	Ø 25 mm	49 mm	
		Rod③	Ø 30 mm	43 mm	
		Rod④	Ø 34 mm	35 mm	
		Joint①	Ø 37,5 mm	82,3 mm	
		Joint②	Ø 42,8 mm	82,6 mm	
		Joint③	Ø 47,1 mm	82,3 mm	
3-9798-03	5355-0250	Handle	Ø 51,5 mm	115,2 mm	0.493
		Rod①	Ø 20 mm	126 mm	
		Rod②	Ø 25 mm	120 mm	
		Joint①	Ø 25 mm	37,5 mm	

3-9798-04	5355-0280	Rod①	Ø 20 mm	90 mm	0.8655
		Rod②	Ø 25 mm	84 mm	
		Rod③	Ø 30 mm	78 mm	
		Rod④	Ø 34 mm	70mm	
		Joint①	Ø 37,5 mm	82,3 mm	
		Joint②	Ø 42,8 mm	82,6 mm	
		Joint③	Ø 47,1 mm	82,3 mm	
		Handle	Ø 51,5 mm	115,2 mm	
3-9798-05	5355-0300	Rod①	Ø 20 mm	114 mm	0.741
		Rod②	Ø 25 mm	108 mm	
		Rod③	Ø 30 mm	100 mm	
		Joint①	Ø 37,5 mm	82,3 mm	
		Joint②	Ø 42,8 mm	82,6 mm	
		Handle	Ø 46,9 mm	115,3 mm	
3-9798-06	5355-0450	Rod①	Ø 20 mm	163 mm	1.013
		Rod②	Ø 25 mm	158 mm	
		Rod③	Ø 30 mm	150 mm	
		Joint①	Ø 37,5 mm	82,3 mm	
		Joint②	Ø 42,8 mm	82,6 mm	
		Handle	Ø 46,9 mm	115,3 mm	
3-9798-07	5355-0600	Rod①	Ø 35 mm	137 mm	1.67
		Rod②	Ø 30 mm	151 mm	
		Rod③	Ø 25 mm	151 mm	
		Rod④	Ø 20 mm	161 mm	
		Joint①	Ø 47.1 mm	82,3 mm	
		Joint②	Ø 42,8 mm	82,6 mm	
		Joint③	Ø 37,5 mm	82,3 mm	
		Handle	Ø 46,9 mm	115,3 mm	