





品番	型番	容量	実測値（測定箇所は上図参照）※単位：mm								目盛（mL）		
			a	b	c	d	e	f	g	h	開始	終了	1目盛
4-1624-01	446081	50mL	35	42	38	48	69	35	15	5	4	50	2
4-1624-02	447081	100mL	43.5	53	48	61	85	43	19	6	4	100	2
4-1624-03	440081	250mL	56	70	63	83	120	62	30	10	15	250	5
4-1624-04	441081	500mL	76	92	85	105	140	72	32	12	30	500	10
4-1624-05	442081	1L	93	115	103	137	181	105	45	15	30	1000	10
4-1624-06	443081	2L	121	150	138	170	213	97	47	18	60	2000	20
4-1624-07	444081	3L	135	167	153	188	245	107	50	18	150	3000	50
4-1624-08	445081	5L	165	178	204	227	270	127	50	20	150	5000	50