



Document identification number : V4FV-NIWX-AQEC-DV7S

Certificate of Laboratory Testing

Issue No.	R20200506-0006	Receipt No.	170200022-001	
Test Completion Date	2017-02-14	Date of Receipt	2017-02-01	
Product Name	Dr.Puri Micro-Dust Protection Mask Large (KF94)			
Manufacturing report number		Batch Number	7701	
Product Type, Material, Classification	Health Mask			
Date of Manufacturing	2017-01-31			
Expiration Date	2020-01-30 3Years from date of manufacture			
Client	Company Name	KM Corporation		
	Name	Shin byung soon		
Manufacturer	Address	16-49, Je2gongdan 5-gil, Miyang-myeon, Anseong-si, Gyeonggi-do, Republic of Korea (Tel. No.: Fax No.: E-mail.)		
	Name	KM Corporation		
	Country	Korea, Republic of		
Test Purpose	Etc			
Remarks				
Test Items and Results				
Test Items	Specifications	Test Results	Judgement	Remarks
Appearance	Visual examination	Suitability	Same as above	
Shape test	Visual examination	Suitability	Same as above	
Acid or Alkali / Fluorescent / Pigment	There should be no red color / There should be no fluorescent / There should be no color	Suitability	Suitability	
Formaldehyde	It should not be darker than potassium chromate	Suitability	Suitability	
Tensile strength	More than 10 N	33 N	Suitability	
Particulate filtration efficiency(sodium chloride)	More than 94 %	96.7 / 96.9 / 96.7 / 96.8 / 96.9 / 96.7	Suitability	

* This certificate has been issued on the Internet, you can check whether to forge or modify by the issued number. Furthermore, You can also check it by barcode exploiting document check program for scanner

인력 2020-05-06 16:21:32



문서확인번호(VERIFICATION NO.) : W4PB-R0XE-U7HH-HB2N

Ministry of Food and Drug Safety
Gyeongin Regional Office of Food and Drug Safety
 Building #5, Gwacheon Government Complex, 47 Gwanmoon-ro, Gwacheon-si, Gyeonggi-do, Korea, Tel: 82-2-2110-8000, Fax: 82-2-2110-0801

Certificate of Free Sales

- No. of Certificate : 2020-D1-1712
- └ Exporting (certifying) country : Republic of Korea
- └ Importing (requesting) country : EU

The Ministry of Food and Drug Safety, certifies that the following firm is authorized to manufacture quasi drug under The Pharmaceutical Affairs Act, and the following product(s) is(are) permitted to be freely sold in domestic and overseas markets.

- Applicant (=Product-license holder)
(This certificate shall not be issued to others than the product-license holder)

- Name : Km corporation
- Address : 16-49 Je2gongdan 5-gil, Miyang-myeon, Anseong-si, Gyeonggi-do, Republic of Korea
- Registered No : 1020

Name of product	Ingredient(s)	No. and date of product-license, comments
Dr.Puri Micro-Dust Protection Mask KF94 (Large, Medium)	Mask	19,2017-04-28

* Attached, if necessary : approved product information (O)

Issued date : APR. 29, 2020 (Certificate No.2020-D1-1712)

Certified by Moon Byung Ho

Moon Byung Ho

Director General Services Division
 Gyeongin Regional Food & Drug Administration

- 1 -



You can verify the Certificate through VERIFICATION NO. on the webpage(<https://nedrug.mfds.go.kr/pbp/CCBSC03/certificate>) or by checking the barcode with the mobile scanner App (MaSmartDetector).