

# お客様のご用途に合った最適なシャフトジェネレーターをお選びください。

処理範囲 リットル	10 ml				50 ml				100 ml				500 ml				減圧下での使用	周速度 (m/s)	懸濁液の最終 粒度 (μm)	乳液の最終粒 度 (μm)
	10 ml	50 ml	100 ml	500 ml	2 l	10 l	20 l	50 l	2 l	10 l	20 l	50 l								
A	S 10 N-5 G	■												No	6.0	5-25	1-10			
	S 10 N-8 G	■	■											No	9.6	5-25	1-10			
	S 10 N-10 G	■	■	■										No	11.9	5-25	1-10			
	S 10 D-7 G-KS-65	■	■											No	7.5	10-50	5-20			
	S 10 D-7 G-KS-110	■	■											No	7.5	10-50	5-20			
B	S 18 N-10 G	■	■	■										No	9.8	10-50	1-10			
	S 18 N-19 G	■	■	■	■									No	16.6	10-50	1-10			
	S 18 D-10 G-KS	■	■	■										No	8.8	10-50	5-20			
	S 18 D-14 G-KS	■	■	■	■									No	12.4	10-50	5-20			
	S 25 N-8 G	■	■											No	8	10-50	1-10			
C	S 25 N-10 G	■	■	■										No	9.8	10-50	1-10			
	S 25 N-18 G	■	■	■	■									No	16.6	10-50	1-10			
	S 25 KV-18 G	■	■	■	■									Yes	16.6	10-50	1-10			
	S 25 NK-19 G	■	■	■	■									No	16.6	10-50	1-10			
	S 25 N-25 G	■	■	■	■	■								No	22.3	15-50	1-10			
D	S 25 KV-25 G	■	■	■	■	■								Yes	22.3	15-50	1-10			
	S 25 N-25 F	■	■	■	■	■								No	23.6	5-25	1-5			
	S 25 KV-25 F	■	■	■	■	■								Yes	23.6	5-25	1-5			
	S 25 D-10 G-KS	■	■	■										No	8.8	10-50	5-20			
	S 25 D-14 G-KS	■	■	■	■									No	12.4	10-50	5-20			
E	S 50 N-G 45 G	■	■	■	■	■	■							No	18.8	40-100	10-30			
	S 50 N-G 45 M	■	■	■	■	■	■	■						No	21.2	25-50	5-20			
	S 50 N-G 45 F	■	■	■	■	■	■	■	■					No	20.9	10-30	1-10			
S 65 KG-HH-G 65 G	■	■	■	■	■	■	■	■	■				Yes	21.9 (28.8: T 65 digital)	25-75	5-25				
S 65 KG-HH-G 65 M	■	■	■	■	■	■	■	■	■	■			Yes	21.9 (28.8: T 65 digital)	25-50	5-15				
S 65 KG-HH-G 65 F	■	■	■	■	■	■	■	■	■	■	■		Yes	21.9 (28.8: T 65 digital)	5-20	1-10				



T 10  
basic

A



T 18  
digital

B



T 25  
digital

C



T 50  
digital

D



T 65  
basic

E



T 65  
digital