


# チューブクイックジョイント 6-670-01 ~ 04

## 寸法図



## 寸法表

(単位:mm)

コード	-01	-02	-03	-04
型番	S	M	L	LL
A	57.8	65.5	70.4	75.4
B	34.3	38.2	42.5	42.8
C	34.3	38.1	40.3	42.7
D	3.1	4.0	6.0	7.8
E	4.0	5.6	7.8	10.1
F	5.1	6.6	8.8	11.2
G	6.1	7.5	10.0	12.4
H	6.0	8.1	10.6	13.3
I	15.7	15.6	17.5	18.3
J	7.5	7.5	9.5	9.5
K	6.4	6.2	8.1	8.2
L	11.0	10.9	12.5	12.5
M	15.7	15.6	17.5	18.2
N(内径)	1.8	3.0	4.0	5.5