



Specification

1.14291.4000 Acetonitrile isocratic grade for liquid chromatography LiChrosolv®

Specification		
Purity (GC)	≥ 99.8	%
Identity (IR)	conforms	
Evaporation residue	≤ 4.0	mg/l
Water	≤ 0.05	%
Acidity	≤ 0.0005	meq/g
Alkalinity	≤ 0.0002	meq/g
Transmission (at 195 nm)	≥ 70	%
Transmission (at 200 nm)	≥ 90	%
Transmission (from 240 nm)	≥ 98	%
Filtered by 0.2 µm filter		

Jeannette David

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.