




designed for scientists



## R 3004.1 Blade stirrer

/// Data Sheet

Generates a radial flow. For drawing the material to be mixed from the top and bottom. For aerating of low viscosity media and for homogenizing solids in liquids. Used at medium to high speeds.



designed for scientists

### Technical Data

Stirrer diameter [mm]	50
Shaft diameter [mm]	10
Shaft length [mm]	565
Material in contact with medium	AISI 316L / 1.4571
Weight [kg]	0.383
Permissible ambient temperature [°C]	5 - 40

