

1. Part Number: 8015A-170S

2. Photograph:



3. Product Information:

| | |
|---|---|
| Description: | SteriWare Micro-Spatula - 170mm - Sterile |
| Use: | Single use, disposable powder spatula |
| Material of Construction: | High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> • FDA acceptable CFR 177.1640 • EU Regulations 10/2011 compliant • EC Regulations 1935/2004 compliant • BSE/TSE free |
| Method of Construction: | Single piece, injection moulded |
| Spatula Length: | 175mm |
| Spatula Blade Length: | 40mm |
| Blade Width: | 13mm |
| Nominal Weight: | 5g |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Number of Samplers per Box: | 100 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |
| ATEX: | Suitable for use in ATEX rated areas – Contact Sampling Systems for a copy of the report to confirm suitability for your application |

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Jonathan Fenner

Position: Director

Date: 21-May-19