











FLAME TIP Compatibility List

Pipette brand manufacture		Eppendoer Research Plus (Single)			Eppendoer Research Plus (Multi)			Eppendorf Explorer(Multi)	Gilson Piepetman (Single)					Gilson Piepetman (Multi)		Sartorius (Single)					Sartorius (Multi)			Thermo Finpipette (Single)					
Model / range		0.5-10µl	20-200µl	100-1000µl	0.5-10µl	10-100µl	30-300µl	120-1200µl	50-1200l	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl	P10, 0.5-10µl	P300, 20-300µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	10-100µl	30-300µl	20-200µl	100-1000µl	
Volume	Product No.	FilterTip																											
 0.1-10µl	FTS252010	●			●					●	●					●		●							●				
 0.5-10µl long	FTS233010	●			●					●	●					●		●							●				
 2-20µl	FTS252020		●			●	●					●	●	●		●		●	●	●	●				●		●	●	
 10-100µl	FTS252100		●			●	●					●	●			●			●	●	●				●	●	●	●	●
 10-200µl	FTS231200		●			●	●					●	●			●				●	●	●			●	●	●	●	●
 10-300µl	FTS233300		●			●	●					●	●			●					●	●			●	●	●	●	●
 100-1000µl	FTS252000			●				●							●									●				●	●
 100-1000µl long	FTS371000			●				●							●									●				●	●