

**KPL-U10065-JP-GCQ (WHITE)**
KPL-U20065-JP-GCQ (BLACK)**SPECIAL FEATURES**

- Washable medical keyboard
- 1.8 meter long USB cable with waterproof USB cap
- 3 standard keyboard LEDs (Num, Scroll, Caps Lock)
- Can be cleaned with common detergents

SPECIFICATION

Number of keys:	102(JP-Layout)	Cable length:	1.8 meter USB cable
Keyswitch technology:	Membrane	Operating temperature:	0°C ~ +55°C
Switching cycles:	Approx. 2 mill. (per key)	Storage temperature:	-15°C ~ +60°C
Housing material:	Plastic	Dimensions:	441(L) x 143(W) x 35.3(T) mm
Housing color:	White	Weight:	Approx. 590g (cable included)
Interfaces:	USB	Certification:	CE, RoHS, IEC 60601-1-2:2014

ORDER INFORMATION

Cat. No.	Color	Layout
KPL-U10065-JP-GCQ	White	JP
KPL-U20065-JP-GCQ	Black	JP

Other layouts, configurations and interfaces on request.

SUPPLEMENT PRODUCTS

Cat. No.	Product description	Color
MSI-U10030	Waterproof medical touch scroll mouse	White
MSI-U20030	Waterproof medical touch scroll mouse	Black




GETT.ASIA

GETT Asia Ltd.

WASHABLE KEYBOARDS

JP LAYOUT



CLEANING INSTRUCTION

1) USE DISINFECTANT SPRAYS

Thoroughly spray all surface area and between keys. Wipeaway excess moisture prior to use. It is recommend to use the disinfectant sprays in line with the manufacturer's specifications.

2) USE DISINFECTANT WIPES

Use a disinfectant wipes in-line with manufacturer's specification to clean keyboard. Exert slight pressure and wipe the keyboard surface with circular movement.

For further details please refer to below link:

<https://gett-group.com/en/gett-service/cleaning-and-material-resistances>

GETT Gerätetechnik GmbH

Mittlerer Ring 1
08233 Treuen
Germany

Tel: +49 37468 660 0
E-Mail: info@gett-group.com
Internet: www.gett-group.com

GETT Asia Ltd.

Room 806-7, Kowloon Plaza
No. 485 Castle Peak Road
Kowloon, Hong Kong

Tel: +852 3461 9789
E-Mail: info@gett-asia.com
Internet: www.gett-asia.com

GETT Assembly (Asia) Co., Ltd.

4F & 5F, Building 1, No.97, Fuyong Area 107 Road
Baishixia Community, Fuyong Street
Baoan District, Shenzhen, China

Tel: +86 755 2323 6358
E-Mail: info@gett-asia.com
Internet: www.gett-asia.com