


ナフロン®プラグ付分液ロート（玉はパイレックス® ガラス製） 7-188-01～08

外形寸法図



< 寸法表 >

型番	-01	-02	-03	-04	-05	-06	-07	-08
容量	50ml	100ml	200ml	250ml	300ml	500ml	1000ml	2000ml
A	53mm	66mm	81mm	83mm	89mm	106mm	132mm	165mm
B	238mm	254mm	293mm	305mm	317mm	354mm	385mm	440mm
C	$\frac{15}{\text{§}} / 25$	$\frac{15}{\text{§}} / 25$	$\frac{24}{\text{§}} / 30$	$\frac{24}{\text{§}} / 30$	$\frac{24}{\text{§}} / 30$	$\frac{24}{\text{§}} / 30$	$\frac{29}{\text{§}} / 32$	$\frac{29}{\text{§}} / 32$
D	13.5mm	13.5mm	16mm	16mm	16mm	18mm	18mm	24.5mm
E	120mm	120mm	130mm	130mm	140mm	150mm	155mm	150mm
F	8mm	8mm	10mm	10mm	10mm	12mm	12mm	16mm