


外形寸法図



< 寸法表 >

型番	-09	-10	-11	-12	-13	-14	-15
容量	50ml	100ml	200ml	250ml	300ml	500ml	1000ml
A	55mm	60mm	73mm	75mm	83mm	98mm	126mm
B	204mm	245mm	292mm	292mm	305mm	350mm	415mm
C	§ ¹⁵ / ₂₅	§ ¹⁵ / ₂₅	§ ²⁴ / ₃₀	§ ²⁴ / ₃₀	§ ²⁴ / ₃₀	§ ²⁴ / ₃₀	§ ²⁹ / ₃₂
D	13.5mm	13.5mm	16mm	16mm	16mm	18mm	18mm
E	60mm	60mm	70mm	70mm	70mm	80mm	90mm
F	8mm	8mm	10mm	10mm	10mm	12mm	12mm