


Screw-Cap Tubes

The Screw-Cap Tube is commonly used in sample degasification. The vacuum quality that it can maintain is $>10^{-4}$ torr. For better vacuum, please check our Pressure/Vacuum Tube and Quick Pressure Valve Tube. Each Screw-Cap Tube comes with one PTFE/Silicone Septum.

Note:

PTFE/ Rubber septums are inert to most solvents and many corrosive materials but not recommended for multiple punctures.

PTFE/Silicone septums are inert to most organic solvents and compounds but not recommended for strongly corrosive materials. The septum remains reliable after multiple punctures.



Screw-Cap Sample Tube

Catalog No.	MHz Rating	O.D.	Length
335-TR-7	600	3mm	7"
335-TR-8	600	3mm	8"
328-TR-7	500	3mm	7"
328-TR-8	500	3mm	8"
328-TR-9	500	3mm	9"
307-TR-7	300	3mm	7"
535-TR-7	600	5mm	7"
535-TR-8	600	5mm	8"
535-TR-9	600	5mm	9"
528-TR-7	500	5mm	7"
528-TR-8	500	5mm	8"
528-TR-9	500	5mm	9"
507-TR-7	300	5mm	7"
507-TR-8	300	5mm	8"
507-TR-9	300	5mm	9"
513-7TRA-7	500	10mm	7"
513-7TRA-8	500	10mm	8"
513-7TRA-9	500	10mm	9"

Screw-Cap Sample Tube Replacement Parts

Catalog No.	Description	Fits Tube with O.D.	Package Qty.
TR-LR-01	PTFE/rubber septum	4 and 5mm	36
TR-LR-05	PTFE/rubber septum	10mm	36
TR-LR-07	PTFE/rubber septum	12, 13, 15, 16, and 18mm	36
TR-LS-01	PTFE/silicone septum	4 and 5mm	36
TR-LS-03	PTFE/silicone septum	7.5 and 8mm	36
TR-LS-05	PTFE/silicone septum	10mm	36
TR-LS-07	PTFE/silicone septum ²	12, 13, 15, 16, and 18mm	36
TR-SC-01	Solid Cap	4 and 5mm	12
TR-SC-05	Solid Cap	10mm	12
TR-SC-07	Solid Cap	12, 13, 15, 16, and 18mm	12
TR-SC-09	Solid Cap	20mm	12
TR-OC-01	Open Cap	4 and 5mm	12
TR-OC-03	Open Cap	7.5 and 8mm	12
TR-OC-05	Open Cap	10mm	12
TR-OC-07	Open Cap	12, 13, 15, 16, and 18mm	12