

## ProCulture® Spinner Flasks

- Impeller eliminates dead spots and extends above media surface to increase aeration
- Sidearm openings large enough for 50 mL pipettes
- Flask base has central bump to prevent material from accumulating at bottom center
- Complete assembly



VWR Cat. No.	Capacity	Center Neck Size	Side Arm Thread
76335-212	125ml	70mm	GL32
76335-214	250ml	70mm	GL32
76335-216	500ml	100mm	GL45
76335-218	1 L	100mm	GL45
76335-220	3 L	100mm	GL45
76335-222	6 L	100mm	GL45
76335-224	8 L	100mm	GL45
76335-226	15 L	100mm	GL45
76335-228	36 L	100mm	GL45

## ProCulture® Shaker Flasks

### Delong Shaker Flasks (A)

- Long neck reduces risk of splashing that could reach the closure

### Beaded Rim Shaker Flasks (B)

- Beaded neck rim ensures greater resistance to chipping

### Fernbach Shaker Flasks (C)

- Ideal for cultures that require large surface area to volume ratio



VWR Cat. No.	Capacity	O.D.	Height	Neck O.D.	Baffle Style
<b>ProCulture® Delong Shaker Flasks (A)</b>					
76323-662	250ml	81mm	138mm	38mm	No Baffles
76323-664	500ml	101mm	181mm	38mm	No Baffles
76323-666	1000ml	129mm	285mm	38mm	No Baffles
76323-668	2000ml	160mm	285mm	38mm	No Baffles
76323-674	125ml	66mm	122mm	25mm	Side Baffles
76323-676	250ml	81mm	138mm	38mm	Side Baffles
76323-678	500ml	101mm	181mm	38mm	Side Baffles
76323-680	1000ml	129mm	285mm	38mm	Side Baffles
76323-682	2000ml	160mm	285mm	38mm	Side Baffles
76323-684	4000ml	206mm	305mm	38mm	Side Baffles
76323-686	125ml	66mm	122mm	25mm	Deep Side Baffles
76323-688	250ml	81mm	138mm	38mm	Deep Side Baffles
76323-690	500ml	101mm	181mm	38mm	Deep Side Baffles
76323-692	1000ml	129mm	285mm	38mm	Deep Side Baffles
76323-694	2000ml	160mm	285mm	38mm	Deep Side Baffles
76323-696	4000ml	206mm	305mm	38mm	Deep Side Baffles

VWR Cat. No.	Capacity	O.D.	Height	Stopper Size	Neck O.D.	Baffle Style
<b>ProCulture® Beaded Rim Shaker Flasks (B)</b>						
76323-698	125ml	66mm	122mm	#5	25mm	Side Baffles
76323-700	250ml	81mm	138mm	#6	38mm	Side Baffles
76323-702	500ml	101mm	181mm	#7	38mm	Side Baffles
76323-704	1000ml	129mm	285mm	#8.5	38mm	Side Baffles

VWR Cat. No.	Capacity	O.D.	Height	Neck Size	Stopper Size	Baffle Style
<b>ProCulture® Fernbach Shaker Flasks (C)</b>						
76323-670	2800ml	200mm	230mm	Beaded Rim	#13	No Baffles
76323-672	2800ml	200mm	230mm	Beaded Rim	#13	Side Baffles
76323-586	2800ml	200mm	230mm	GL45 threaded	GL45*	Side Baffles
76323-588	2800ml	200mm	260mm	Delong	—	Bottom Center Baffles

\*GL45 cap included