




System contents and spare parts for reaction units 100-500 ml, DN 60

Cat. No.: 6 2601 37, Cat. No.: 6 2601 49, Cat. No.: 6 2601 58

Pos.	System contents (pc.)	Description	Cat. No.
1a	1	Reaction vessel 100 ml, DN 60	6 0301 37
1b	1	Reaction vessel 250 ml, DN 60	6 0301 49
1c	1	Reaction vessel 500 ml, DN 60	6 0301 58
2	1	Lid, center neck NS 29/32, 2 side necks NS 14/23	6 1031 41
3	1	Stirrer guide, cone NS 29/32	6 3060 04
4	1	Glass stirrer, shaft with lateral PTFE paddle	6 3070 01
5	1	Dropping funnel, 100 ml, NS 14/23, with graduation	4 3737 11
6	1	Dimroth condenser, jacket length:160 mm, NS 14/23	5 2402 01
7	1	Adapter 90°, cone NS 14/23, with olive	5 0950 01
8	1	Hollow stopper, pointed bottom, NS 14/ 23	1 4000 14
9	4	Forked clamps NS 14/23, stainless steel	1 1035 14
9a	1	Forked clamp NS 29/32, stainless steel	1 1035 29
10	1	Tension brace DN 60	6 2465 06
11	1	Support, complete with connector	6 8000 00
12	1	Sealing ring, DN 60, silicone	6 2480 06
13	2	Plastic olives, PP, with screw cap GL 14	1 3212 14
14	1	Spare valve plug, 10 mm, PTFE, black	6 0613 10