

## System contents and spare parts for one-way reaction/distillation units 6, 10, 20 l, DN 150

Cat. No.: 6 2654 85, Cat.No.: 6 2654 88 and Cat. No.: 6 2654 94



Pos.	System contents (pc.)	Description	Cat. No.
1a	1	Reaction vessel 6 l, DN 150	6 0314 85
1b	1	Reaction vessel 10 l, DN 150	6 0314 88
1c	1	Reaction vessel 20 l, DN 150	6 0314 94
2	1	Flat flange lid, DN 150, 1 center neck, 3 side necks	6 1044 65
3	1	Stirrer guide, cone NS 45	6 3062 06
4a	1	Double-propeller stirrer for 6 l	6 3310 06
4b	1	Double-propeller stirrer for 10 l	6 3310 10
4c	1	Double-propeller stirrer for 20 l	6 3310 20
5	1	Hinged coupling, metal	6 4005 01
6	1	Motor, 230 V	6 9000 01
7a	1	Dosing funnel 1000 ml, NS 29/32 (for 6 and 10 l)	6 3500 70
7b	1	Dosing funnel 4000 ml, NS 29/32 (for 20 l)	6 3500 79
8	1	Vent tube	6 3505 04
9	1	Bend 135°, 2 NS 29/32	6 3550 04
10	1	Distilling head 120°	6 3551 14
11	1	Temperature probe Pt 100, 85 mm, NS14/23	6 5000 05
12a	2	Flange support, stainless steel, KF 15 (for 6 and 10 l)	6 2485 15
12b	2	Flange support, stainless steel, KF 25 (for 20 l)	6 2485 25
13a	2	PTFE-seals KF 15 (for 6 and 10 l)	6 2488 15
13b	2	PTFE-seals KF 25 (for 20 l)	6 2488 25
14a	2	Tube adapter KF 15, thread M 16 x 1 (for 6 and 10 l)	5 1081 10
14b	2	Tube adapter KF 25, thread M 16 x 1 (for 20 l)	5 1082 16
15	1	Coil condenser	5 2550 29
16	1	Bend, NS 29/32, S 29	6 3550 10
17	1	Intermediate receiver 0.5 l	6 5150 05
18	1	Final receiver 2.0 l	6 5150 10
19	2	Ventilation valve 0-2.5 mm, NS 29/32	6 5200 00
20	1	Connection valve 10 mm	6 5200 05
21	10	Forked clamps NS 29/32	1 1035 29
22	1	Forked clamp S 29	1 2329 15
23	1	Support DN 150, complete with connector	6 2478 01
24	1	O-ring, silicone, DN 150	6 2480 15
25a	1	Valve plug, 10 mm, PTFE, for 6 l	6 0613 10
25b	1	Valve plug, 20 mm, PTFE, for 10 and 20 l	6 0612 20