



LABORATORY

LABORATORY AND MEDICAL BRUSHES

Dun bristles - nylon - goat hair

Our laboratory brushes are designed to be used for the manual cleaning of hoses, test tubes, spouts, measuring cylinders, pipettes and nozzles.

On request, we can make brushes in any dimensions, diameters and with any bristle materials (see the Material information section above).

Brush end

Wire shaft with eyelet



TEST TUBE BRUSH

Image 1



TEST TUBE BRUSH WITH WOOL TIP

Image 2





PIPETTE BRUSH

Image 3



Image	Application	Article number	Dimensions in mm			
			A	B	C	D
1	Test tube brush	50300401	270	100	120	30
1	Test tube brush	50300403	270	100	120	30
2	Test tube brush	50300301	270	100	130	30
2	Test tube brush	16500301	270	80	105	20
2	Test tube brush	16500303	270	80	105	20
2	Test tube brush	70150301	280	70	85	15
2	Test tube brush	70150303	280	70	85	15
2	Test tube brush	70100301	230	50	60	10
3	Pipette brush	25000904	400	100	-	5x10x5
3	Pipette brush	52800104	490	190	-	4x20x4
3	Pipette brush	10780104	450	180	-	25/12
3	Pipette brush	10(øø)0101	500	100	-	3-20
3	Pipette brush	10(øø)0103	500	100	-	3-25
3	Pipette brush	10(øø)0104	500	100	-	3-20
4	Spout or hose brush	50260101	270	90	-	26
4	Spout or hose brush	50260103	270	90	-	26
4	Spout or hose brush	16500101	270	80	-	20
4	Spout or hose brush	16500103	270	80	-	20
4	Spout or hose brush	70100101	230	50	-	10
4	Spout or hose brush	70100103	230	50	-	10
4	Spout or hose brush	70150101	270	70	-	15
4	Spout or hose brush	70150103	270	70	-	15
5	Nozzle brush	19100101	140	45	-	10
5	Nozzle brush	19100103	140	45	-	10
5	Nozzle brush	19100104	140	45	-	10
6	Pipette brush	26000404	250	20	25	20x10



Twist wire	Bristles	Brush tip	Brush end	Image
Galvanised	Dun bristles	Tip bent	Wire with eyelet	1
Galvanised	Nylon	Tip bent	Wire with eyelet	1
Galvanised	Dun bristles	With wool tip	Wire with eyelet	2
Galvanised	Dun bristles	With wool tip	Wire with eyelet	2
Galvanised	Nylon	With wool tip	Wire with eyelet	2
Galvanised	Dun bristles	With wool tip	Wire with eyelet	2
Galvanised	Nylon	With wool tip	Wire with eyelet	2
Galvanised	Dun bristles	With wool tip	Wire with eyelet	2
Galvanised	Goat hair	Rhombus-shaped centre	Wire with eyelet	3
Galvanised	Goat hair	Rhombus-shaped	Wire with eyelet	3
Galvanised	Mane goat hair	Two-part	Wire with eyelet	3
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	3
Galvanised	Nylon	Slightly angled	Wire with eyelet	3
Galvanised	Goat hair	Slightly angled	Wire with eyelet	3
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Dun bristles	Slightly angled	Wire with eyelet	4
Galvanised	Nylon	Slightly angled	Wire with eyelet	4
Galvanised	Goat hair	Slightly angled	Wire with eyelet	5
Galvanised	Goat hair	Slightly angled	Wire with eyelet	5
-	Goat hair	Conical, angled	Wire with eyelet	6



Dimensions A B C D

Should you wish us to make a special design to your requirements, please specify: overall length, bristle length, diameter, bristle thickness and bristle material

Please note that we use only high quality materials. If you are unable to find a brush appropriate for your intended application in our catalogue, contact us; we can design a suitable brush for you and produce it at a competitive price.

SPOUT OR HOSE BRUSH

Image 4



NOZZLE BRUSH

Image 5



PIPETTE BRUSH

Image 6



26000404