

BagPipet® & BagTips®

Fixed-volume pipet and disposable straws



- **Fully autoclavable**
- **Economical**
- **Straw ejector**
- **Easy to calibrate**
- **No contamination**
- **No blocked straws**



BagPipet® is an automatic pipet which intakes preset sample volumes of 0.1, 0.9 or 1 mL, depending on the volume button chosen. It can make serial dilutions in tubes, 1 mL in 9 mL, or dispense 1 mL or 0.1mL of sample in a Petri dish.

BagPipet® is made of light metal and is fully autoclavable.

Fitted with a straw ejector, BagPipet® uses sterile **BagTips®** straws in packets of 25, which can be re-closed using a **BagClip®** closing clip. The gamma-sterilized BagTips® are highly rigid and low cost.

- 2 models of BagTips®:**
- **REGULAR (19 cm / 7.5")**.
 - **JUMBO (24 cm / 9.5")**.



BagPipet® is particularly well suited for intaking thick liquid or food samples after blending in a BagMixer® with a BagFilter® or BagPage® filter bag.

Please send me: a quotation for **BagPipet® & BagTips®**.
 information on **BagSystem®**.
 free **BagFilter®** sample.

Name

Company.....

Address.....

.....

Zip code / City.....

Country.....

Tel..... Fax.....

e-mail..... web.....



International : 30. Ch Bois Arpents, 78860 St NOM FRANCE
 Tel. +33(0)1 34 62 62 61 / Fax +33(0)1 34 62 43 03
 info@interscience.fr / www.interscience.fr

USA/Canada: 199Weymouth St. ROCKLAND, MA 02370 USA
 Tel. +1 781 792 2133 / Fax +1 781 792 2134
 info@intersciencelab.com / www.intersciencelab.com

www.interscience.fr