

# MAGLab Series™

## Create lab storage space — off your bench top

With lab bench top space at a premium, common lab products such as pipettors, test tubes, Kimwipes®\* and miscellaneous accessories can be easily organized with Heathrow Scientific's magnetic lab storage system, the MagLab. MagLab includes four magnetic pieces: "MagRack" that holds (9) 15ml and (4) 50ml tubes; "Mag Wipe" that holds a standard size (11 cm x 21cm) Kimwipes; box; "MagPette" that holds two pipettors (universal design fits most common brands); and "MagBox" that holds miscellaneous items such as pens/pencils, notes, etc. Each rack is made of robust lightweight ABS plastic in bright and fun translucent colors and features four powerful magnets for a secure hold onto most metal surfaces or create your own holding point with the metal plate included for each segment. The MagLab Series is available individually or as a set of four.

Part No.	Description	L x W x H		Color
HS23500	maglab series, 4/pk	20.5 x 15.4 x 4 in	52.2 x 39.3 x 10.2 mm	■ ■ ■ ■ *
HS23501	magrack	5.1 x 4.4 x 2.5 in	13.1 x 11.2 x 6.4 mm	Teal ■
HS23502	magwipe	5.1 x 3.5 x 3.8 in	13.0 x 9.1 x 9.8 mm	Green ■
HS23503	magpette	5.1 x 4.4 x 2.9 in	13.1 x 11.2 x 7.6 mm	Purple ■
HS23504	magbox	5.1 x 3 x 4 in	13.0 x 7.8 x 10.2 mm	Orange ■

\*Assorted colors: Teal, Green, Purple and Orange

Kimwipes® is a registered trademark of Kimberly-Clark Corp.



Includes self-adhesive metal plates that adhere to any surface

Heathrow Scientific®